# **B.Sc. CHEMISTRY LAB MANUAL**

3rd Semester

Prepared By Pure & Applied Science Dept. Chemistry

# MIDNAPORE CITY COLLEGE

# MIDNAPORE CITY COLLEGE

# Department of Pure and Applied Sciences B.Sc. in Chemistry Semester: III Paper Title: Physical Chemistry-I (Practical) Paper Code: CEMHMJ03P

#### **Course contents:**

1. Calibration of various apparatus such as burette, pipette, volumetric flux, measuring cylinder etc.

2. Concepts of molecular weight, equivalent weight, various concentration terms, primary and secondary standard solutions and their preparations with proper explanations of types of apparatus to be used to prepare these solutions.

Prepare the following solutions

(i) Primary standard oxalic acid solution by accurate weighing

(ii) Primary standard potassium dichromate solutions by accurate weighing

(iii) Secondary standard sodium hydroxide solutions by approximate weighing

(iv) Secondary standard sodium thiosulphate solutions by approximate weighing

(v) Secondary standard potassium permanganate solution by approximate weighing

(vi) Secondary standard sulphuric acid, hydrochlroric acid and acetic acid solutions of various strength from the supplied concentrated solutions.

3. Study of kinetics of acid-catalyzed hydrolysis of an ester

4. Study of kinetics of decomposition of H2O2

5. Determination of specific rotation of cane sugar using polarimeter

6. Study the kinetics of inversion of cane sugar using polarimeter.

#### MIDNAPORE CITY COLLEGE

2

# Experiment 1: Calibration of various apparatus such as burette, pipette, volumetric flask, measuring cylinder, etc.

# **Introduction:**

You have studied the use of various volume measuring devices such as measuring cylinders, burettes, pipettes, and volumetric flasks previously. Measuring cylinders are used for approximate measurement of volume for various purposes and do not require any calibration. Burettes, pipettes, and volumetric flasks are used for accurate quantitative measurements. When using apparatus of standard quality. it is not desirable to attempt calibration of such apparatus for three reasons. Firstly. the process of calibration is time-consuming. Secondly, due to the possibility of introduction of errors in calibrating the apparatus, the apparatus thus calibrated by an inexperienced beginner may be less accurate than the apparatus calibrated by the manufacturer. Thirdly, the error permitted in a determination may take care of any inaccuracies in the calibration of the apparatus. However, if the apparatus being used is not of standard quality. it is advisable to check the calibration of the apparat& before performing any measurements with such an apparatus. In this experiment. you will perform the calibration of a pipette, a burette, and a volumetric flask.

#### **Objectives:**

After performing this experiment, you should be able to:

• Justify the need for calibration of apparatus

• State the steps involved in the calibration of volumetric apparatus, calibrate burettes, and volumetric flasks.

#### Principle

Calibration of a burette and pipette is carried out by measuring the mass of a known volume of water delivered by the burette or the pipette. By knowing the temperature and the density of water at that temperature, the true volume delivered by the burette or the pipette is calculated. This will be the volume delivered by the burette or the pipette and not the one which is etched on the apparatus. The volumetric flasks are calibrated by filling them using a calibrated pipette.

#### Requirements

Apparatus: Analytical balance -1, Beaker -1, Burette -1, Conical flask- 4, Cork/stoppers- 4, Pipette -1, Thermometer -1, Volumetric flask -1

3

Chemicals: Distilled water

#### Procedure

Calibration of a Pipette:

Take a clean 100 cm3 (mL) conical flask fitted with a cork or stopper. Weigh the conical flask with cork accurately on an analytical balance and record its mass in your notebook. Pipette out 10/20/25 cm3 (equal to the capacity of the pipette u have) of distilled water from a beaker into the conical flask and weigh the flask with its contents. Record the temperature of the water with a thermometer. Calculate the mass of water delivered from the difference between the two readings. Repeat the procedure 3-4 times. Take the average mass and divide it by the density of water at room temperature to calculate the actual capacity of the pipette.

Temperature	Density	Temperature	Density
°C	g/ cm <sup>3</sup>	°C	g∕ cm³
15	0.99913	26	0.99682
16	0.99897	27	0.99655
17	0.99880	28	0.99627
18	0.99862	29	0.99598
19	0.99843	30	0.99568
20	0.99823	31	0.99537
21	0.99802	32	0.99506
22	0.99780	33	0.99473
23	0.99757	34	0.99440
24	0.99733	35	0.99406
25	0.99708		

Values of the density of water at different temperatures are below

#### **Calibration of a Burette**

Take a clean 100 cm3 (mL) conical flask fitted with a cork. Measure its mass accurately to the nearest milligram on an analytical balance and record the mass. Clean the burette thoroughly and grease the stop-cock properly. Fill up the burette with distilled water and note the room temperature. Ad-just the water level to the zero mark. Touch the tip of the burette to the wall of a beaker to remove any adhering drop. Take the initial reading of the water level in the burette. Run 10 cm3 of distilled water into the conical flask. Touch the tip of the burette I to the wall of the flask. Cork the flask and weigh it accurately on an analytical balance and record the mass. The difference between this mass and the initial mass of the flask gives the mass of

water delivered. Calculate the volume of water actually delivered by dividing the mass by the density of water. Repeat the process with 20 cm3, 30 cm3, 40 cm3, and 50 cm3 water, thus calibrating the burette over its entire length.

#### **Calibration of a Volumetric Flask**

Volumetric flasks may also be calibrated like that described for burettes and pipettes. But their large volume requires the use of an oversized analytical balance. Therefore, the volumetric flasks are calibrated by the procedure described below.

Take a clean 250 cm3 volumetric flask. Using a calibrated 25 cm3 pipette, carefully transfer 25 cm3 of distilled water 10 times in succession, to the flask. If the level of water in the flask does not correspond with the mark on the neck of the flask. paste a label with its upper edge on a level with the bottom of the meniscus of water. If you want to find the true volume of the water in the flask up to the graduation mark shown on the flask, you may do so with the help of a calibrated burette by transferring water exactly up to the graduation mark.

#### **Observations**

#### **Calibration of a Pipette**

- Mass of the empty corked conical flask = W<sub>1</sub>g
- ii) Mass of the corked conical flask + distilled water =  $W_2g$
- iii) Temperature of water  $\dots = t^{o}C$
- iv) Density of water at  $t^{\circ}C = d$

#### **Calibration of a Burette**

- Temperature of water = t°C
- ii) Density of water at  $t^{o}C = d$
- iii) Mass of the empty corked conical flask =  $w_1 g$
- iv) Mass of the corked conical flask + 10 cm<sup>3</sup> distilled water =  $w_2g$
- v) Mass of the corked conical flask + 20 cm<sup>3</sup> distilled water =  $w_3 g$
- vi) Mass of the corked conical flask + 30 cm<sup>3</sup> distilled water =  $w_4g$
- vii) Mass of the corked conical flask + 40 cm<sup>3</sup> distilled water =  $w_5 g$
- viii) Mass of the corked conical flask + 50 cm<sup>3</sup> distilled water =  $w_6 g$

#### **CALCULATIONS:**

#### Calibration of pipette

Mass of water delivered by the pipette =  $(W_2 - W_1)g$ 

Volume of water delivered by the pipette =  $\frac{W_2 - W_1}{d} = X \text{ cm}^3$ Hence, the pipette is calibrated to deliver =  $X \text{ cm}^3$ 

#### **Calibration of burette**

Mass of 10 cm<sup>3</sup> water delivered by the burette =  $(w_2 - w_1) g$ 

Volume of 10 cm<sup>3</sup> water delivered by the burette =  $\frac{w_2 - w_1}{d} = x_1 \text{ cm}^3$ 

Error in volume

 $=(x_1-10)$  cm<sup>3</sup>

 $= \frac{(x_1 - 10) \times 100}{x_1}$ 

% Error

Similarly calculate the calibrated volume and error in burette for  $20 \text{ cm}^3$ ,  $30 \text{ cm}^3$ ,  $40 \text{ cm}^3$  and  $50 \text{ cm}^3$  and tabulate the results as given below.

S.No.	Volume of Water delivered	Volume of water as found by calibration	Error in cm <sup>3</sup>	% Error
1	10 cm <sup>3</sup>	x1 cm3	$(x_1 - 10)$	$100(x_1 - 10)/x_1$
2	20 cm <sup>3</sup>	x <sub>2</sub> cm <sup>3</sup>	$(x_2 - 20)$	$100(x_2 - 20)/x_2$
3	30 cm <sup>3</sup>	x <sub>3</sub> cm <sup>3</sup>	$(x_3 - 30)$	$100(x_1 - 30)/x_1$
4	40 cm <sup>3</sup>	x4 cm3	$(x_4 - 40)$	$100(x_4 - 40)/x_4$
5	50 cm <sup>3</sup>	x <sub>5</sub> cm <sup>3</sup>	$(x_{5} - 50)$	$100(x_5 - 50)/x_5$

#### Summary:

In this experiment, you would have performed the following: calibrated a pipette, calibrated a burette, and calibrated a volumetric flask.

Experiment 2: Concepts of molecular weight, equivalent weight, various concentration terms, primary and secondary standard solutions and their preparations with proper explanations of types of apparatus to be used to prepare these solutions.

Prepare the following solutions

(i) Primary standard oxalic acid solution by accurate weighing

(ii) Primary standard potassium dichromate solutions by accurate weighing

(iii) Secondary standard sodium hydroxide solutions by approximate weighing

(iv) Secondary standard sodium thiosulphate solutions by approximate weighing

(v) Secondary standard potassium permanganate solution by approximate weighing

(vi) Secondary standard sulphuric acid, hydrochlroric acid and acetic acid solutions of various strength from the supplied concentrated solutions.

#### **Experiment 3: Study of kinetics of acid-catalyzed hydrolysis of ester**

#### **Theory:**

The acid catalysed hydrolysis reaction of ester can be represented as

 $R_1COOR_2 + H_2O + H^+ \rightleftharpoons R_1COOH + R_2OH$ 

The rate equation is represented as

$$\frac{d[R_1 COOR_2]}{dt} = k[R_1 COOR_2][H^+][H_2O] = k_1[R_1 COOR_2],$$

where  $k_1 = k[H^+][H_2O] = constant$ , (in large excess of water, [H<sub>2</sub>O] remains constant and H<sup>+</sup> being the catalyst [H<sup>+</sup>] is constant)

Therefore the reaction becomes first order w.r.t ester. Integration of this equation with the boundary conditions, when t= 0,  $[R_1COOR_2] = C_0$  and at time t,  $[R_1COOR_2] = C$ , gives

$$K_1 = (2.303/t) \log (C_0/C)$$

The progress of the catalysed reaction may be studied by withdrawing measured volume of aliquot from the reaction mixture at different intervals of time and titrating with standard alkali solution using phenolphthalein indicator. The volume of alkali required at any instant is equivalent to be the sum of weak acid (produced as a result of hydrolysis) and the acid used as catalyst (a constant quantity).

If  $V_0$ ,  $V_n$  and  $V_{\infty}$  be the volumes of alkali required for the same volume of aliquot at the beginning, at time t= t<sub>n</sub> and at the end of the reaction (at infinite time) respectively then,

$$C_0 \propto (V_{\infty} - V_0)$$
 and  $C \propto (V_{\infty} - V_n)$   
Then,  $k_1 = \frac{2.303}{t_n} \log \frac{(V_{\infty} - V_0)}{(V_{\infty} - V_n)}$ 

To avoid the measurement of  $V_0$  the equation may be represented in the form

$$\Delta t_n = t_n - t_1 = \frac{2.303}{k_1} \log \frac{(V\omega - V1)}{(V\omega - Vn)}$$

Where  $V_n$  and  $V_1$  are the volumes of alkali required at times  $t_n$  and  $t_1$  respectively.

Plot of  $log \frac{(V\omega - V1)}{(V\omega - Vn)}$  vs  $\Delta t_n$  will give a straight line passing through the origin and the slope  $k_1$  may be determined by

 $K_1$ =slope× 2.303

# **Apparatus required:**

- 1) 100 mL dry conical flask-1
- 2) 500 mL bottle for NaOH.
- 3) 2 mL pipette; 5 mL pipette-1 each
- 4) 250 mL conical flask-5
- 5) Sufficient ice cold water; water bath
- 6) Stop watch

#### Chemicals required: Methyl acetate, NaOH, and Phenolphthalein

#### **Procedure:**

- 1) Prepare 250 mL of approximately 0.1(N) NaOH solution.
- 2) Prepare 100 mL 1(N) HCl solution
- 3) Pipette out 50 mL of prepared HCl solution (catalyst solution) in a 100 mL dry conical flask, and add 5 mL of ester using a pipette. Start the stopwatch at the time of half discharge. Mix the solution thoroughly by a swirling motion.
- 4) At 5-7 minute intervals, take a 2 mL aliquot and add to 50 mL ice-cold water taken in a 250 mL conical flask. Note the time of half discharge. Titrate rapidly against the prepared~ 0.1(N) NaOH solution taken in a burette, using phenolphthalein as an indicator. Take at least 6 readings.
- 5) The remaining solution is heated at about 60°C in a water bath~for 40 minutes with an air condenser fitted in the mouth of the conical flask. The solution is allowed to cool to room temperature. Pipette out 2 mL of it in 50 mL of water taken in a 250 mL conical flask and titrate with~0.1(N) NaOH solution, using phenolphthalein as an indicator. The titre value corresponds to V∞.

- 6) Plot a graph of  $log \frac{(V\omega V1)}{(V\omega Vn)}$  vs  $\Delta t_n$  and draw the best fit straight line passing through the origin.
- 7) Calculate the value of ' $k_1$ ' from the slope of the graph.

# **Experimental result:**

1. Room temperature

2. Preparation of 250 mL~0.1(N) NaOH solution: Dissolve 1 g NaOH in 250 mL deionized water.

3. Preparation of 100 mL 1(N) HCl solution: Dissolve~9 mL conc. HCl in 91 mL deionized water.

4. Recording of data for ester hydrolysis:

Time(t)	Time(t <sub>n</sub> )in sec	$\Delta t_n = t_n - t_1$	Volume of NaOH (V <sub>n</sub> mL)	$(V_{\omega}-V_1)$ mL	$(V_{\infty}-V_n)$ mL	$log \frac{(V\omega - V1)}{(V\omega - Vn)}$
			,			

5. Determination of  $V_{\infty}$ :

 $V_{\infty}$  = Volume of NaOH required to standardize 2 mL reaction mixture after heating at about 60°C in a water bath.

Plotting of graph

$\log \frac{(V\omega - V1)}{(V\omega - V1)}$			
$\frac{\log 100}{(V\infty - Vn)}$			
$\Delta t_n(Sec)$			

**Calculation:** Observed rate constant, k<sub>1</sub>=slope×2.303

**Conclusion:** Therefore observed rate constant for the acid catalyzed hydrolysis of methyl acetate at .....<sup>o</sup>C is.....sec<sup>-1</sup>.

#### Experiment 4: Study of kinetics of decomposition of H<sub>2</sub>O<sub>2</sub>

#### **Theory:**

The reaction between H<sub>2</sub>O<sub>2</sub> and KI in dilute H<sub>2</sub>SO<sub>4</sub> medium is represented as

$$H_2O_2 + 2KI + H_2SO_4 = I_2 + 2H_2O + K_2SO_4$$

Suggested reaction mechanism is

$$H_2O_2+I^{-}=H_2O+OI^{-}$$
 ....(slow)  
 $OI^{-}+2H^{+}+I^{-}=H_2O+I_2$  .....(fast)

The first step being the slowest step is the rate determining step. The differential rate becomes

 $\frac{d[H_{2O_2}]}{dt} = k[H_2O_2][I]$ , where k is the rate constant and other terms have their usual significance.

Therefore, kinetically the reaction is of second order (overall). But if the concentration of iodide ion is kept constant [by adding Na<sub>2</sub>S<sub>2</sub>O<sub>3</sub> continuously whereby the S<sub>2</sub>O<sub>3</sub><sup>-</sup> ions react with the liberated iodine (I<sub>2</sub>) and regenerate I<sup>-</sup> according to the equation, I<sub>2</sub>+2 S<sub>2</sub>O<sub>3</sub><sup>2-</sup> = S<sub>4</sub>O<sub>6</sub><sup>2-</sup>+2I<sup>-</sup>] the reaction becomes kinetically first order w.r.t.[ H<sub>2</sub>O<sub>2</sub>] only. Under this condition the rate equation may be expressed as

$$\frac{d[H_{2O_2}]}{dt} = k_1[H_2O_2], \text{ where } k_1 = k[I^-]$$

If the initial (at time t=0) concentration of  $H_2O_2$  is a and at time t,  $[H_2O_2] = (a-x)$ , (where, x= amount of  $H_2O_2$  reacted), integration of the above equation gives

$$K_1 = \frac{2.303}{t} \log \frac{a}{a - x}$$

If both 'a' and (a-x) are represented in terms of their  $Na_2S_2O_3$  equivalent, we have  $V_0 \propto a$  (where  $V_0$  is the titre value of thiosulphate for the total iodine liberated by a fixed volume of  $H_2O_2$  solution of concentration 'a' after complete decomposition).

And  $V_t \propto x$  (where  $V_t$  is the titre value of the same thiosulphate solution for the iodine liberated by the same H<sub>2</sub>O<sub>2</sub> solution undergoing reaction at time t).

The rate equation may then be expressed as

$$K_{1} = \frac{2.303}{t} \log V_{0}/V_{0}-V_{t}$$
  
Or, log V\_{0}/V\_{0}-V\_{t}=k\_{1}t/2.303

A plot of  $V_0/V_0-V_t$  vs. t will give a straight line passing through the origin and from the slope  $k_1$  may be determined,  $k_1$ = slope× 2.303.

MIDNAPORE CITY COLLEGE

### **Apparatus required:**

- 1) 500 mL conical flask-2
- 2) 500 mL glass bottle for thiosulphate-1
- 3) Burette-1
- 4) 250 mL volumetric flask for KI-1
- 5) 10 mL pipette
- 6) Watch glass
- 7) Stop watch

**Chemicals required:** H<sub>2</sub>O<sub>2</sub>, Sodium thiosulphate, KI, 12N H<sub>2</sub>SO<sub>4</sub>, 1% Ammonium molybdate solution, Starch

#### **Procedure:**

1. Prepare 100 mL '2 volume' H<sub>2</sub>O<sub>2</sub> solution. in a glass bottle.

2. Prepare 250 mL of 0.1(N) sodium thiosulphate solution in a glass bottle. Fill the burette with sodium thiosulphate solution.

3. Take 10 mL 12 (N)  $H_2SO_4$  in a 500 mL conical flask. Add 10 mL of prepared  $H_2O_2$  solution followed by addition of 50 mL of water. Then add approximately 2 gm soli KI. Add 2 mL 1% ammonium molybdate solution. Cover the conical flask with a watch glass and keep it in dark for 1 minute. Add about 50 mL of water and titrate the liberated iodine with sodium thiosulphate solution using starch as an indicator. Record the volume of thiosulphate (V<sub>0</sub>).

4. Take 250 mL of 0.4% KI solution (dissolve 1 gm KI in 250 mL of deionized water) in a 500 mL conical flask. Add 15 mL 12(N) sulphuric acid and 5 L of freshly prepared starch solution. Add 10 mL of prepared  $H_2O_2$  solutions and at the time of half discharge start stop watch. Run immediately sodium thiosulphate solution from the burette from the burette into the mixture in excess to discharge the blue colour. Wait for appearance of blue colour and record the time (t) of reappearance of colour, volume of thiosulphate used (V<sub>t</sub>). Again run thiosulphate into the mixture in excess to discharge the blue colour, record time (t) of reappearance of colour, and total volume of thiosulphate used (V<sub>t</sub>). Take at least 6 readings.

5. Plot log  $[V_0/(V_0-V_t)]$  vs. t in millimeter graph paper and draw the best fit straight line passing through the origin. Calculate 'k<sub>1</sub>' value from the slope of the graph.

#### **Experimental result:**

- 1) Room temperature
- Preparation of 250 mL~0.1(N) sodium thiosulphate solution: Dissolve~6.2 gm Na<sub>2</sub>S<sub>2</sub>O<sub>3</sub> in 50 mL deionized water.
- 3) Determination of V<sub>0</sub>:

 $V_0$  = titre value of thiosulphate for the total iodine liberated by 10 mL of '2 volume'  $H_2O_2$  solution after complete decomposition.

4) Titre value of thiosulphate at different times:

Time(t)	Time (t) in sec	Volume of	$(V_0-V_t) mL$	log
		$S_2O_3^{2-}$ used (V <sub>t</sub>		$V_0/V_0-V_t$
		mL)		

#### 5. Plotting of graph:

$\log V_0/V_0-V_t$			
t(sec)			

**Calculation:** Observed rate constant,  $k_1$ =slope× 2.303

**Conclusion:** Therefore, the observed rate constant for the decomposition of  $H_2O_2$  at .....<sup>9</sup>C is ......Sec<sup>-1</sup>.

# Experiment 5: Determination of specific rotation of cane sugar using polarimeter

#### **Object:**

To determine the specific rotation of cane sugar solution using Polarimeter.

#### **Specific Rotation:**

Specific rotation is the rotation produced by a column in a solution of an optically active substance one centimeter long, the concentration of the substance being one gram per c.c.

#### **Apparatus Used:**

Polarimeter, Electric bulb with housing, sugar, balance, graduated cylinder and thermometer.

#### Formula:

The specific rotation S of the plane of polarization of sugar dissolved in water.

~ /1

	$S = \theta / lxc = \theta x v / l x m$
Where,	c = Concentration of the solution m/v =gm/cc
	$\theta$ = Rotation produced in degree
	l = Length of the tube in cm
	m = Mass of sugar in g dissolved in water
	v = Volume of sugar solution

#### Method:

1. Weigh 5 gm of sugar and carefully dissolve it is water to make up to 100cc of solution. The solution should be well mixed by pouring it from one jar to another so that it becomes of one uniform concentration. This gives a 5% solution of sugar.

2. Clean the glass tube, fill it with water and close the end, see that there are no air-bubbles in the tube and no water is sticking on the outside of the end glasses. Put the tube in position and rotate the analyzer. A till the field of view is completely and uniformly dark i.e. both halves are equally dark. Read the position of analyzer on the circular scale, provided for the purpose. A further rotation through 1800 will again produce an extinction.

3. Remove the tube, empty it, rinse it with solution and then fill it completely, removing all bubbled. Clean the ends and put the tube in position. The two halves of the field will appear unequally dark.

4. Rotate the analyzer to produce extinction again and read the scale. This difference gives the rotation produced. Get a second position of extinction (equal darkness) by a further rotation of 1800 and take mean for the rotation produced.

5. Now, to 50 c.c. of solution add 50 c.c of water and mix thoroughly. This gives a 5% solution and mix thoroughly. This gives a 5% solution. Take readings with it and then mix 50 c.c. of this solution with 50 c.c. of water again. This gives a 2.5% solution. Take readings with this also. Repeat with 1.25% solution.

13

#### Table for **\theta**

Value of one division of mains scale = No. of division on Vernier scale = Least count of Vernier analyzer =

#### Table -1

S.NO.	% solution i.e. reading of analyzer with water						
	Clockwise			A	$a = \frac{X+Y}{X+Y}$		
	M.S.	V.S.	Total 'X'	M.S.	V.S.	Total 'Y'	2
1.							
2.							
3.							

#### Table – 2

S.NO.	Concentration of solution	% solution i.e. reading of analyzer with water				$Mean \\ b = \frac{X' + Y'}{2}$	$(\boldsymbol{\theta}) = \boldsymbol{a} - \boldsymbol{b}$		
		(	Clockwise Anticlockwise						
		M.S.	V.S.	Total	M.S.	V.S.	Total		
				<b>'X'</b>			<b>'Y'</b>		
1.	5% solution								
2.	2.5% solution								
3.	1.25% solution								

### **Calculations:**

1. Draw a graph between  $\theta$  and concentrations. The graph will be a straight line as shown in Fig.1.

2. Find out the value of  $\theta$  for particular value of concentrations.

3. Calculate the value of specific rotations S as

 $S = \theta x v/l x m ...^{o}/cm/gm/cc.$ 





#### **Result:**

- 1. The specific rotation of sugar solution (Solvent water) =  $\dots$  o/cm/gm/cc at °C at  $\lambda$  = 5880 Ao.
- 2. Standard value of specific rotation of water =  $66.5^{\circ}$ .
- 3. Standard value of specific rotation of Glucose (Solvent water) =  $52^{\circ}$ .

#### **Precautions and Source of Error:**

- 1. The polarimeter tube should be well cleaned.
- 2. Water and sugar used should be dust free.

3. Whenever a solution is changed, risen the tube will new solution under concentration or clean the tube will water and dry it.

- 4. There is should be no bubble inside the tube.
- 5. The position of analyzer should be set accurately.

6. Readings should be taken when halves of the field of view becomes equally illuminated.

# Experiment 6: Study the kinetics of inversion of cane sugar using polarimeter

#### **Theory:**

if a compound when placed in the path of plane of polarization it is said to be optically active and this property of rotating the plane of polarized light is called optical activity. If a compound rotates the plane of polarization in clock clockwise direction, facing the beam, it is called dextrorotatory, and if it rotates in an anti-clockwise direction, it is called laevorotatory. A solution of sucrose is optically active and dextrorotatory, while glucose and fructose are dextrorotatory and laevorotatory.

#### **Principle:**

An important characteristic of reaction is that sucrose is dextrorotatory, whereas glucose and fructose are dextrorotatory and laevorotatory. Further, the laevorotation of fructose is more (being-92°) than the dextro rotation of glucose (being -52.50), so that the mixture as a whole is laevorotatory. It means the product has a resultant levorotatory effect. Thus on hydrolysis, the dextrorotatory sucrose gradually changes into the levorotatory mixture. It is for this reason that the reaction is called the Inversion of Sucrose. The inversion of cane sugar is a pseudo-first-order reaction and represented by

#### $C_{12}H_{22}O_{11} + H_2O {\rightarrow} C_6H_{12}O_6 + C_6H_{12}O_6$

And the angle of rotations of Sucrose is (+66.6°), while that of Glucose and Sucrose are (-52.5°) and (-92°) respectively.

#### Materials:

Polarimeter, sample tube, flasks, beakers, sodium lamp and stirrer.

#### **Chemicals:**

Sugar solution and HCl.

#### **Procedure:**

Firstly, prepared 5% sucrose solution and 0.1 N HCl solution in separate flasks and placed them in water bath to maintain the temperature. Flasks should be closed with stopper to avoid any fumes. Then mixed equimolar solution of both HCl and sucrose. Temperature should be maintained. Took sample tube, cleaned its lens and poured sample solution in it form its one end and placed the sample tube in polarimeter, applied light source from one end and note the reading from the main scale and Vernier scale at time 0,7.18, 27.05 min and infinity, and noted the specific rotations at different time intervals.

#### **Observation and Calculation:**

Reading of polarimeter at any time  $= r_t$ 

Reading of polarimeter at zero time =  $r_0$ =+24.09

Reading of polarimeter at infinite time =  $r_{\infty} = -10.74$ 

Sr no.	Time (min)	Rotations (degree)	$\frac{t^{p} - r\infty}{rt - r\infty}$	$\mathbf{K} = \frac{2.303}{t} \log \frac{v^{*} - r\infty}{r t - r\infty}$ (min <sup>-1</sup> )
1	7.18	+21.41	1.083	$1.1 \times 10^{2}$
2	18	+17.74	1.223	1.1×10 <sup>2</sup>
3	27.05	+15.0	1.353	1.1×10 <sup>2</sup>

#### Table: Record of data to determine the rate constant for inversion of cane sugar.

The average value of rate constant  $k = 1.1 \times 2^{-2} \text{ min}^{-1}$ .

#### **Result:**

As the value of k is constant, it indicates that the inversion of cane sugar is first order reaction.

# Paper Title: Organic Chemistry-III (Practical)

# Paper Code: CEMHMJ04P

#### **Course contents:**

#### **Experiment: Qualitative Analysis of Single Solid Organic Compounds**

- a) Detection of special elements (N, S, Cl, Br) by Lassaigne's test
- b) Solubility and classification (solvents: H2O, 5% HCl, 5% NaOH and 5% NaHCO3)

c) Detection of the following functional groups by systematic chemical tests: aromatic amino (-NH2), aromatic nitro (-NO2), amido (-CONH2, including imide), phenolic –OH, carboxylic acid (-COOH), carbonyl (-CHO and >C=O).

d) Melting point of the given compound

# **Qualitative Analysis of Single Solid Organic Compounds**

# (A) Physical Characteristics

- (i) State of aggregation: solid
- (ii) Texture: Crystalline or amorphous
- (iii) Colour:
- (iv) Odour:

#### (B) Melting Point:

#### (C) Solubility Test

Cold H <sub>2</sub> O	Hot H <sub>2</sub> O	dil. HCl	dil. NaHCO <sub>3</sub>	dil. NaOH

#### **Conclusion:**

- From solubility and melting point inference can be drawn for the given organic sample.
- If the sample is soluble in dil. NaOH. The compound may contain the functional group
   COOH or phenolic –OH or –COOH and phenolic –OH both.
- As the sample effervesces with NaHCO<sub>3</sub> solution i.e the compound is soluble in NaHCO<sub>3</sub> solution then the sample may contains –COOH group with or without other functional groups.

17

- As the sample is soluble in only dil. HCl then the compound is amine with or without other functional groups.
- As sample is soluble in both dil. HCl and dil. NaOH the compound may be aminophenol or amino acid.
- As the sample is not soluble in dil acid or alkali but only soluble in concentrated H₂SO₄.
   The compound may contain −NO₂ and or >C=O group.

# (D) Detection of Special Elements

### Lassaigne's Test:

Take a small piece of dry sodium metal in a fusion tube and heat it gently till the metal melts or fuses. Add equal quantity of compound to this fused metal. Heat it gently then strongly till it becomes red hot. Plung the red hot tube in 10 ml distilled water taken in a porcelain dish, crush the fusion tube completely. Carry out more fusion in the similar way. Boil the extract for five minutes, reduce the volume to about 5 ml and filter. Perform following test using this filtrate.

Test	Observation	Inference
Test for Nitrogen: Filtrate + freshly	Prussian blue or green ppt or	N-present
prepared FeSO4 soln. + boiled, cooled +	colour.	
acidified with dil. H2SO4 soln.		
<b>Test for Sulphur:</b> 1 ml of filtrate + 1 drop	A violet or purple coloration	S-present.
of sodium nitroprusside solution.		
<b>Test for Halogen:</b> 1 ml of filtrate + conc.	White ppt. which is soluble in	Cl-present.
HNO <sub>3</sub> , boiled, cooled, + AgNO <sub>3</sub> soln.	dil. NH4OH but reappears	
	when acidified with HNO <sub>3</sub>	

**Conclusion:** From the above experiment, is concluded the supplied organic sample contain N, S, and Cl as special elements.

# (E) Systematic analysis of the functional groups:

#### 1. Detection of Nitrogenous functional groups

Experiment Observation Inference
----------------------------------

Test for aromatic-NH2 group: Dye test:	Red or orange-red dye	Ar-NH <sub>2</sub> group
Solution-I: Sample is dissolved in dil. HCl and		confirmed.
cooled in an ice-bath.		
Solution-II: Dil. NaNO2 soln. is cooled in an ice-		
bath.		
Solution-III: Alkaline $\beta$ -napthol soln. is cooled		
in an ice-bath.		
The mixture of solution-I & II is cooled in an		
ice-bath. After some time alkaline $\beta$ -napthol		
soln. is added to the mixture of solution-I & II.		
Test for -NO <sub>2</sub> group: Mulliken and Barker's	Silver mirror or black or	Ar-NO <sub>2</sub> group
test:	grey ppt.	confirmed.
Alcoholic soln. of sample + solid NH4Cl + Zn-		
dust boiled for few minutes, cooled allowed to		
stand for 5 minutes, and filtered. Filtrate +		
Tollen's reagent, warmed in a water bath.		
Test for -CONH <sub>2</sub> group: Hydrolysis test:	Characteristic smell of	-CONH <sub>2</sub> group
sample + Conc. NaOH solution, heated.	NH3 which turns	confirmed.
	mercurous nitrate paper	

# 2. Detection of Non-Nitrogenous functional groups:

Experiment	Observation	Inference
Test for -COOH group: (i) NaHCO3 Test:	Effervescence of	-COOH group may be
Saturated solution of sodium bicarbonate under heating, cool + pinch of sample.	CO <sub>2</sub> gas.	present.
(ii) Esterification test: 1.0 g sample + 5 ml	A sweet smell of	-COOH group is
ethanol + 2-3 concentrated $H_2SO_4$ , boiled in a water bath for 10 minutes. The mixture is then poured into 150 ml water.	ester.	present and confirmed.
Test for phenolic-OH group: (i) FeCl <sub>3</sub> Test:	Violet or green	Phenolic-OH group
alcoholic or aqueous solution of sample + FeCl <sub>3</sub>	or blue	may be present.
solution.	colouration.	

(ii) Back dye test: Solution-I: small amount of	Orange-red dye.	Phenolic-OH group
aniline is dissolved in dil. HCl, cooled in ice-		present and
bath.		confirmed.
Solution-II: NaNO <sub>2</sub> solution, cooled in ice-bath.		
Solution-III: alkaline solution of sample, cooled		
in ice-bath.		
The solution-I & II are mixed and kept for some		
time. The solution-III is added to the mixture of		
solution-I & II.		
Test for carbonyl group: 2 g sample is	Yellow or	Carbonyl (>C=O or -
dissolved in minimum quantity of ethanol + 5 cc	orange ppt.	CHO) group is present
saturated solution of 2,4-		and confirmed.
dinitrophenylhydrazone, warmed, scratched		
with glass rod.		
Test for -CHO group: (i) Fehlings's test:	Yellow or red	-CHO group present
sample + Fehlings's solution (Fehlings-I &	ppt.	and confirmed.
Fehlings-II), heated in water bath.		
(ii) Tollen's Test: sample + Tollen's reagent,	Silver mirror or	-CHO group present
heated in water bath.	grey ppt.	and confirmed.

#### **Highlights:**

1. During Lassaigne's test a dark coloured sodium extract is obtained either due to incomplete fusion of the organic compound or when it is taken in excess of the sodium taken. Again if sodium is not used in excess during the experiment, on pouring the fusion tube into water no sodium hydroxide will be formed and the solution will not be alkaline. So green precipitate of ferrous hydroxide will not be observed during Prussian blue test and the test for N will fail.

2. If the solution is not alkaline sodium nitroprusside test will also fail.

Above two tests can be performed by adding sodium hydroxide solution to the filtrate if the solution is not alkaline.

3. The solution will give a black precipitate on addition of ferrous sulphate if sulphur is present in the organic sample.

20

4. The ferrous sulphate solution should be freshly prepared. A little amount of ferrous sulphate is washed with distilled water to remove ferric sulphate is more soluble than ferrous sulphate. Next the sulphate is dissolved to get freshly prepared ferrous sulphate solution. A dilute FeSO<sub>4</sub> on standing is oxidized to basic Ferric Sulphate which becomes useless in analysis.

5. The addition of Ferric chloride solution and HCL is avoided during Prussian blue test as it forms a green precipitate or colour. Concentrated sulphuric acid, in turn, will give a fine blue precipitate or colouration as concentrated sulphuric acid will oxidise the excess ferrous sulphate present in solution to ferric sulphate and those ferric icons will yield Prussian blue in fine blue colour.

6. It is essential to boil the sodium extract with HNO<sub>3</sub> to remove CN- as voilatole HCN. Otherwise it will give give ppt of AgCN similar to AgCl and misguide the analysis.

 $NaCN + AgNO_3 = AgCN + NaNO_3$ 

 $NaCN + HNO_3 = HCN + NaNO_3$ 

# Paper Title: Preparation of Household Chemicals (Practical) Paper Code: CEMSEC03

#### **Course contents:**

#### Part A:

- (i) Preparation of liquid disinfectant
- (ii) Preparation of Liquid Detergent
- (iii) Preparation of Dish wash Liquid Soap
- (iv) Preparation of Hand Sanitizer
- (v) Preparation of simple organophosphates, phosphonates and thiophosphates

Part B: Field visit and submission of the report

#### (i) Preparation of liquid disinfectant

#### HOW TO PREPARE DISINFECTANT SOLUTIONS

- Use the concentration of sodium hypochlorite available in the country, as indicated on the label of the container, to mix the concentration of disinfection solution according to the intended use, following the indications in the table below.
- Use protective equipment such as gloves, mask, eye protection, and apron for mixing the solutions.
- Prepare solutions daily. If you don't need a liter per day, make adjustments to the amounts in the table below.
- Use a dedicated container for the prepared solutions, and mark (label) the container with the concentration clearly.
- Do not reuse packaging of edible products to prepare the disinfectant solution to avoid accidental poisoning.
- Add the quantity of sodium hypochlorite to the quantity of water as indicated in the table below.
- > Do not mix sodium hypochlorite with soaps or cleaners.

MIDNAPORE CITY COLLEGE

Store sodium hypochlorite out of the reach and sight of children, in a safe place that is not exposed to light and heat.

Sodium Hypochlorite					Disinfec	tant Soluti	on Con	icentratio	'n			
concentration		0.059	61	0.1%				0.5%				
(as per label)	For ha	or hand washing and laundry floors, and cleaning utensils and mortuary				For spills of body fluids, (blood, vomiting, etc)						
	So hypo	dium chlorite	Qua W	ntity of ater	Soo Нуро	dium chlorite	Qua W	ntity of /ater	So Hypo	dium chlorite	Qua W	ntity of /ater
1%	50 ml	1/5 cup	1 L	1 qt	100 ml	2/5 cup	11	1 qt	500 ml	1 pt	1 L	1 qt
3%	16 ml	3 tsp	11	1 qt	30 ml		11	1 qt	300 ml	½ pt	11	1 qt
4%	13 ml	1/2 oz	11	1 qt	25 ml	1 02	11	1 qt	250 ml	1 cup	1L	1 qt
5%	10 ml	2 tsp	11	1 qt	20 ml	4 tsp	1L	1 qt	100 ml	2/5 cup	1L	1 qt
10%	Sml	1 tsp	11	1 qt	10 ml	2 tsp	11	1 qt	50 ml	1/5 cup	11	1 qt

➤ Use sodium hypochlorite products without added fabric softeners, colors, or perfumes.

# What is Dettol and what is it used for?

Dettol is a brand of multi-purpose liquid disinfectant that is mostly used domestically in our various homes or places of work, for diverse purposes. The liquid disinfectant is mostly used for bathing and washing, it can also be used when cleaning utensils, cleaning wound areas of the skin, and lots more.

Dettol helps to kill harmful germs and other bacteria, or as the case may be, and helps to prevent the spread of germs and harmful bacteria, and protect against illness or infections that can be caused by germs.

Apart from usage in homes, Dettol Liquid Disinfectant can be used in various schools, hospitals, and even churches, etc. to prevent the spread of germs and other diseasecausing bacteria. Apart from bathing and washing, use it too when scrubbing floor tiles, it works wonders too!

Description of Chemicals used in the production of Dettol Liquid Disinfectant, and their respective functions.

# 1. Texapon:

The Taxapon is a jellylike and pasty substance, it functions as a foaming agent in the production. The Texapon is a detergent that produces micelles that allows non-polar substance such as oils to be dissolved in water.

#### 2. Phenol:

Phenol, also known as hydroxyl benzene is a colorless aromatic compound with a very sweet smell (it sometimes, comes in brownish color too). It is one of the disinfectants in used in the production of Dettol; however, in small concentrations. Phenol can also be used as a mouth wash.

#### .3. Chloroxylenol

This is also known as para-chloro-meta-xylenol (or PCMX), it is the main antiseptic agent in the production of Dettol, and helps to kill germs and for disinfection of the skin. Also used together with alcohol for cleaning and disinfecting surgical instruments.

In higher concentrations, Chloroxylenol may cause some form of irritation to the skin and inflammation, so care must be taken to ensure the measurement during the production of the Dettol.

#### 4. Isopropyl Alcohol:

Isopropyl Alcohol also known as IPA, is also one of the vital ingredients used in the production of Dettol Liquid Disinfectant. It is a colorless compound, with a very strong smell. Its powerful antiseptic effects make Dettol have the ability to fight and kill bacteria, fungi, and viruses, especially in wound treatment.

#### 5. Pine oil:

Pine oil is an extract from Pine trees. It is also an antiseptic and disinfectant substance with a very pleasant smell and can also be used as a healing oil for cold. Pine oil also helps to relieve skin inflammation. When in water, it emulsifies (i.e. forms a cloudy solution.) It makes the Dettol form that cloudy solution when dropped in water.

The use of Pine oil is very verse and not limited to the above-mentioned, it can also be used for bathing antiseptic soaps, bath oils, etc. It can be effective too in the healing of arthritis.

# Required Measurement For The Production Of 5 Liters Of Dettol Liquid Disinfectant

- 1. TEXAPON 1/4 kg
- 2. PINE OIL 1/2 liter
- 3. PHENOL 1/8 liter
- 4. CHLOROXYLENOL 1/4 liters
- 5. ISOPROPYL ALCOHOL (IPA) 2 liters
- 6. WATER (Distilled) 2litres
- 7. COLORANT as desired.

Note, emplace of water, you may use methanol or ethanol; or the mixture of both (i.e. water plus methanol). We only need them to properly dissolve other constituents in the production of Dettol Liquid Disinfectant

More so, if you desire to produce more volumes of the Dettol Liquid Disinfectants, all you need do is to play around with your ratios. For example, you can double the above measurements for a 10 Liters production of the Dettol Liquid Disinfectant, and so on. Consider what we have as the basic measurement.

# Production Procedure For The Dettol Liquid Disinfectant.

# Safety Tips During Production:

- Please do note that the production should never be so urgent or so important that you don't take the time to carry it out safely.
- Make sure you put on the right protections before commencing the production. The need for use of rubber hand gloves, nose masks, face shields, etc., cannot be overemphasized.
- Carry out all production and an open and well-ventilated environment. Also, make sure it is well illuminated.
- In the event of a chemical splash to the skin or other part of the body exposed, quickly wash with clean running water, then proceed immediately to see medical personnel for further checks.
- · Please do not taste any chemical, or directly inhale it.
- Please keep all chemicals away from the reach of children.

# Equipment & Tools

- Clean Bowl that can contain 5liters volume without spillage. Preferably a plastic bowl or container.
- A wooden or plastic stirrer
- Funnel
- Packaging containers

# Simple Step-By-Step Guide

Step 1: Pour the Texapon into empty clean bowl.

Step 2: Add the Pine oil and stir properly until the Texapon is dissolved completely.

- Step 3: Add Phenol to the mixture and stir properly
- Step 4: Add Chloroxylenol and stir properly. Note the mixture will turn slightly brownish.
- Step 5: Gently pour the amount of water and stir gradually. The mixture will immediately become cloudy or whitish; this is normal.
- Step 6: Now add the Isopropyl Alcohol (IPA). Continue to stir until the mixture turns brownish again, to get a perfect blend.
- Step 7: Dissolve the color (brown) in a little quantity of water and mix well. Then add bit by bit into the main bowl, while stirring continuously until the desired coloration is reached. Note that this should be the last constituent to be added. Stir very well.



MIDNAPORE CITY COLLEGE

# (ii) Preparation of Liquid Detergent

#### Introduction

**Liquid Detergent** is generally a cleaning chemical that comes in the form of liquid. It is Easy and convenient to use, Liquid Detergents are generally meant for cleaning of clothes. They are one of the most popular and widely used cleaning liquid for Clothes & some other purpose. Liquid Detergent are chemical solutions used as cleaning agents. They are generally formulated using one or more surfactants to improve their cleaning performance. Liquid Detergent have become an essential commodity in our daily lives. Cleaning products play an essential role in our daily lives. By safely and effectively removing soils, germs and other contaminants, they help us to stay healthy, care for our homes and possessions.

#### **Qualities of good Liquid Detergent**

- 1. It must have attractive scent or perfume.
- 2. It must not itch your hand or body.
- 3. Remove germs from the clothes.

#### Formulation

Major Component of liquid detergent is based on surfactants that remove soil from the fabric and clean it. A combination of anionic and non-ionic surfactants is commonly used for optimum performance. In general, anionic are good for particulate soil removal and non-ionic are good for greasy soil removal. Other components of a liquid laundry detergent formula are detailed in this article. Thus, a liquid laundry detergent formula can be as simple and economical as say 5% active surfactant solution in water with preservatives. The formula will foam and do some cleaning of the clothes. From here one can improve this basic formula to any desired level in order to enhance the performance and position it for various categories such as economical, mid-tier, premium, ultra, and 2X, 3X, 4X etc. Liquid Detergent is formulated using batch or continuous process of Liquid making. These cleansing products contain different ingredients that are used to improve their cleaning performance. The surfactant plays an important role in improving the cleansing action of liquid by reducing the surface tension of wash liquid.

#### **Ingredients:**

The important ingredients of Liquid detergent includes- Acid Slurry, caustic soda, AOS, Perfume for fragrance, water.

#### Liquid Detergent Manufacturing Process:

It is a neutralization process of acid slurry by caustic soda keeping PH 7 to 9.

- Firstly transfer water into the stainless steel tank manually or through pump.
- Now add acid slurry into the tank to turn water into desired colour.
- After that it is neutralized by the aqua solution of caustic soda to get the required PH.
- Now AOS is added into the solution, it is used as a foam booster.
- ♦ While urea is used as a buffering agent and also create transparency.
- ✤ Lastly add colours & perfume.
- Filling of Liquid detergent in bottle.
- Packaging of goods.

# **Raw Material Requirement of Liquid Detergent**

- Acid Slurry (alkyl benzene sulphonate (ABS) or linear alkyl benzene sulphonic acid (LABSA)
- Caustic Soda
- ✤ AOS (Alpha Olefin Sulfonate)
- Urea
- Colour
- Perfume
- Corrugated Boxes
- Water

#### **Conclusion:**

After completion of manufacturing process, product is ready to sell in the market. Liquid Detergent are used for daily routine activities. This machine can be installed with low investment & one can earn a good Margin of profit by doing this business.

# (iii) Preparation of Dish Wash Liquid Soap

Dishwashing Liquid is the main agent for removing unwanted germs and bacteria that sticks in plates and other utensils. The experiment's purpose is to produce dishwashing liquid that is affordable yet effective in cleaning, with the process of preparing the 16L of water contained in a beaker. 1000g of SLES was poured in the water and was stirred constantly until it is completely dissolve. 15g Bubble enchancer is then added and stirred to mix it in the solution. After mixing, 15g Foam booster and 15g degreaser is added and mixed to the

solution. Then they dropped 5gcolorant that gave color to the mixture. 15g neutralizer and stabilizer is then added. Sodium Chloride is to adjust the viscosity of the detergent and added scent. After that, it showed that the produced dishwashing liquid is effective and useful as a cleaning agent, it removed the food residues and greases just like the ones sold commercially. The students also observed the active cleansing agents present in the chemical substance used to produce the dishwashing liquid. The only factor that the students suggest is to try other chemical substance that would improve the fragrance and appearance of the dishwashing liquid.

#### INTRODUCTION

Dish washing liquid is one of the most important and the most commonly used household item found in any kitchen. Dishwashing liquid is used primarily for removing food from used dishes and tableware. Heavy soil (large food particles) is generally scraped from the dishes before using. Although its effectiveness appears to be of a chemical nature, it is actually of a physical nature. Dish washing liquid uses the manipulation of surface tension of water which helps break up dirt and grime on dishes. Dishwashing liquid contain active agents called surfactants, that reduce surface tension of water and helps dirt, grime and grease to be separated. This experiment aims to produce a dishwashing liquid made from various chemicals with cleaning properties and to determine the quality of the handmade dishwashing liquid.

# To make your own dishwashing liquid, simply mix the following ingredients in a clean container:

How to make a liquid dishwashing soap?

- ✤ 1 cup of distilled water.
- ✤ 1/4 cup of liquid castile soap.
- ✤ 1 tablespoon of white vinegar.
- ✤ 1 tablespoon of washing soda.
- ✤ 5-10 drops of essential oil (optional, for fragrance)

OR

#### What is the formula for making dishwashing liquid?

In a large bowl, mix together 1 tablespoon (2 g) of Borax, 1 tablespoon (14 g) of washing soda, 2 tablespoons (30 ml) of white distilled vinegar, and ½ cup (120 ml) liquid castile soap. Stir well to ensure that all of the ingredients are well blended.

#### How to make dishwashing liquid at home: Procedure

- > Add your washing soda or baking soda into a bowl with the grated soap. ...
- > Heat water to a boil and pour it into the bowl of washing soda and grated soap.
- > Whisk the mixture until it is blended thoroughly. ...
- > Add a few drops of lemon and lavender essential oil.

#### What makes a good dishwashing liquid?

For clean dishes and a happy planet, a good dishwashing liquid should be hypoallergenic and free of phosphates, and triclosan. It should not be tested on animals and the packaging should be made from recycled materials or come in a bulk container.

#### **INGREDIENTS:**

- Water
- Sodium Hydroxide (Caustic Soda)
- BIO-SOFT® S-118
- STEOL® CS-230
- AMMONYX® LMDO
- Citric Acid, 50% Solution
- Sodium Chloride
- Preservative, Dye, Fragrance, etc.

#### (iv) Preparation of Hand Sanitizer:

This formulation prescribed the requirement for alcohol based hand sanitizers (hand rub) taking into the WHO recommendation. This does not cover non-alcohol based hand sanitizer.

Alcohol-based hand rubs are the only known means for rapidly and effectively inactivating a wide array of potentially harmful microorganisms on hands. WHO recommends alcohol-based hand rubs based on the following factors:

- Evidence-based, intrinsic advantages of fast-acting and broad-spectrum microbicidal activity with a minimal risk of generating resistance to antimicrobial agents;
- Suitability for use in resource-limited or remote areas with lack of accessibility to sinks or other facilities for hand hygiene (including clean water, towels, etc.);
- Capacity to promote improved compliance with hand hygiene by making the process faster and more convenient;
- Economic benefit by reducing annual costs for hand hygiene, representing approximately 1% of extra-costs generated by Health Care-Associated Infection;
- Minimization of risks from adverse events because of increased safety associated with better acceptability and tolerance than other products.

Health-care a setting currently using commercially-available hand rubs should continue to use them, provided that they meet recognized standards for microbicidal efficacy and are well accepted/tolerated by Health Care-Workers. It is obvious that these products should be regarded as acceptable, even if their contents differ from those of the WHO-recommended formulations.

#### Ingredients for 10,000 mL

- Ethanol 96% 8333 mL or
- If Ethanol is used, denaturing is recommended

- Isopropyl Alcohol 99% 7576 mL or
- Isopropyl Alcohol 91% 8242 mL
  - When the concentration of alcohol (e.g., ethanol or isopropyl alcohol) in the starting ingredient is not exact, the calculation should be adjusted accordingly to ensure a final concentration of at least 80% ethanol or 75% isopropyl alcohol
  - Volume of starting ingredient required = (final % alcohol ) × (final volume of preparation) /(starting % alcohol )
  - Methanol contamination of Alcohol
  - Serious safety concern
  - Substantial methanol exposure can result in nausea, vomiting, headache, blurred vision, permanent blindness, seizures, coma, permanent damage to the nervous system or death
  - > Methanol must not be used as an ingredient or as a denaturant
  - Methanol content must not exceed 630 ppm
  - Methanol NMT 630 ppm
  - Benzene NMT 2 ppm
  - Acetaldehyde NMT 50 ppm
  - Acetal (1,1-diethoxyethane) NMT 50 ppm
  - Sum of all other impurities NMT 300 ppm

#### **Composition of Alcohol-Based Hand Rubs Formulations:**

The following two formulations are recommended for local production with a maximum of 50 liters per lot to ensure safety in production and storage.

#### **Formulation I:**

Final concentrations of ethanol 80% v/v, glycerol 1.45% v/v, hydrogen peroxide 0.125% v/v.

Pour the following ingredients as mentioned in table-1 into a graduated flask/ container:

- a. Ethanol (with a purity 96% v/v),
- b. Hydrogen Peroxide 3% v/v,
- c. Glycerol 98% v/v,

Top up the flask to volume required with distilled water or water that has been boiled and cooled; shake the flask gently to mix the content.

#### **Formulation II:**

Final concentrations of isopropyl alcohol 75% v/v, glycerol 1.45% v/v, hydrogen peroxide 0.125% v/v: Pour the following ingredients of required volume as mentioned in table-2 into graduated flask/ container.

a. Isopropyl alcohol (with a purity of 99.8% v/v),

b. Hydrogen Peroxide 3% v/v, c. Glycerol 98% v/v,

Top up the flask to volume required with distilled water or water that has been boiled and cooled; shake the flask gently to mix the content. The % purity of the ingredients may vary. For volume of the ingredients to be taken is calculated by using the following formula.



The examples of calculation by using above formula for volume of 1000 ml, 10000 ml and 50000 ml of preparation of hand sanitizers are given in the table 1 and 2.

Ethanol 96% ml	Glycerol 98% Ml	Hydrogen Peroxide 3% ml	Distilled Water ml	Total Volume Prepared Ml
833	15	42	110	1000
8333	148	417	1102	10000
41667	740	2083	5510	50000

Table -1 Preparation Ethanol Based Hand Sanitizer

Table-2 Preparation of Isopropyl Alcohol based Hand Sanitizer

Isopropyl Alcohol 99.8% ml	Glycerol 98% ml	Hydrogen Peroxide 3% ml	Distilled Water ml	Total Volume Prepared ml
751	15	42	192	1000
7515	148	417	1920	10000
37575	740	2083	9602	50000

Only pharmacopoeia quality reagents should be used.

#### **Production Method**

#### **Volume of Production, Containers**

**10-litre** preparations: glass or plastic bottles with screw threaded stoppers can be used.

**50-litre** preparations: large plastic (preferably polypropylene, translucent enough to see the liquid level) or stainless steel tanks with an 80 to100 liters capacity should be used to allow for mixing without overflowing. The tanks should be calibrated for the ethanol/isopropyl alcohol volumes and for the final volumes of either 10 or 50 liters. It is best to mark plastic tanks on the outside and stainless steel ones on the inside.

#### Preparation

a. The alcohol for the chosen formulation is poured into the large bottle or tank up to the graduated mark.

b. H<sub>2</sub>O<sub>2</sub> is added using the measuring cylinder.

c. Glycerol is added using a measuring cylinder. As the glycerol is very viscous and sticks to the walls of the measuring cylinder, it can be rinsed with some sterile distilled or cold boiled water to be added and then emptied into the bottle/tank.

d. The bottle/tank is then topped up to the corresponding mark of the volume (10-litre or 50-liter) to be prepared with the remainder of the distilled or cold, boiled water.

e. The lid or the screw cap is placed on the bottle/tank immediately after mixing to prevent evaporation.

f. The solution is mixed by gently shaking the recipient where appropriate (small quantities), or by using a wooden, plastic or metallic paddle. Electric mixers should not be used unless "EX" protected because of the danger of explosion.

g. After mixing, the solution is immediately divided into smaller containers (e.g. 1000, 500 or 100 ml plastic bottles). The bottles should be kept in quarantine for 72 hours. This allows time for any spores present in the alcohol or the new or re-used bottles to be eliminated by  $H_2O_2$ .

#### **Quality Control**

If concentrated alcohol is obtained from local production, verify the alcohol concentration and make the necessary adjustments in volume to obtain the final recommended concentration. An alcoholmeter can be used to control the alcohol concentration of the final use solution; H O concentration can be measured by titrimetric (oxidation-reduction reaction by iodine in acidic conditions). A higher level quality control can be performed using gas chromatography and the titrimetric method to control the alcohol and the hydrogen peroxide content, respectively.

MIDNAPORE CITY COLLEGE

Moreover, the absence of microbial contamination (including spores) can be checked by filtration.

#### Labelling

Labels should include the following information:

- Name of producer
- Date of production and batch number

• Composition: ethanol or isopropanol, glycerol and hydrogen peroxide (% v/v can also be indicated) and the following statements:

- For external use only
- Avoid contact with eyes
- Keep out of reach of children

• Use: apply a palmful of alcohol-based hand rub and cover all surfaces of the hands. Rub hands until dry. Flammable: keep away from flame and heat.

#### **Ingredients**

#### a. Hydrogen Peroxide

While alcohol is the active component in the formulations, certain aspects of other components should be respected. All raw materials used should be preferably free of viable bacterial spores. The low concentration of H O is incorporated in the formulations to help eliminate contaminating spores in the bulk solutions and excipients and is not active substance for hand antisepsis. While the use of  $H_2O_2$  adds an important safety aspect, the use of 3-6% of  $H_2O_2$  for the production might be complicated by its corrosive nature and by difficult procurement in some countries. Further investigation is needed to assess H O availability in different countries as well as the possibility of using a stock solution with a lower concentration.

#### **b.** Glycerol

Glycerol is added to the formulation as a humectant to increase the acceptability of the product. Other humectants or emollients may be used for skin care, provided that they are affordable, available locally, miscible (mixable) in water and alcohol, non-toxic, and hypoallergenic. Glycerol has been chosen because it is safe and relatively inexpensive. Lowering the percentage of glycerol may be considered to further reduce stickiness of the hand rubs.

#### c. Other Additives

It is strongly recommended that no ingredients other than those specified here be added to the formulations. In the case of any additions, full justification must be provided together with documented safety of the additive, its compatibility with the other ingredients, and all relevant details should be given on the product label.

A colorant may be incorporated to differentiate the hand rub from other fluids as long as such an additive is safe and compatible with the essential components of the hand rubs. However, the H O in the hand rub may tend to fade any coloring agent used and prior testing is recommended. The

addition of fragrances is not recommended because of the risk of allergic reactions. To further reduce the risk of abuse and to respect cultural and religious sensitivities, product containers may be labelled simply as "antimicrobial hand rub.

#### d. Water

While sterile distilled water is preferred for making the formulations, boiled and cooled tap water may also be used as long as it is free of visible particles.

#### **Production and Storage**

Special requirements apply for the production and stock piling of the formulations, as well as for the storage of the raw materials. Because undiluted ethanol is highly flammable and may ignite at temperatures as low as 10°C, production facilities should directly dilute it to the abovementioned concentration. The flash points of ethanol 80% ( $\nu/\nu$ ) and isopropyl alcohol 75% ( $\nu/\nu$ ) are 17.5°C and 19°C, respectively and special attention should be given to proper storage in tropical climates. Production and storage facilities should be ideally air-conditioned or cool rooms. Open flames and smoking must be strictly prohibited in production and storage areas. Pharmacies and small-scale production centers supplying the WHO-recommended hand rub formulations are advised not to manufacture locally batches of more than 50 liters at a time. For safety reasons, it is advisable to produce smaller volumes and to adhere to local and/or national guidelines and regulations. The production should not be undertaken in central pharmacies lacking specialized air conditioning and ventilation.

#### **Distribution**

To avoid contamination with spore-forming organisms, disposable bottles should preferably be used. To prevent evaporation, containers should have a maximum capacity of 500 ml and possibly fit into a wall dispenser. Leakage-free pocket bottles with a capacity of no more than 100 ml should also be used. The production or re-filling unit should follow norms on how to clean and disinfect the bottles (e.g. autoclaving, boiling, or chemical disinfection with chlorine). Autoclaving is considered the most suitable procedure. Reusable bottles should never be refilled until they have been completely emptied and then cleansed and disinfected.

Chemical disinfection should include soaking the bottles in a solution containing 1000 ppm of chlorine for a minimum of 15 minutes and then rinsing with sterile/cooled boiled water. After thermal or chemical disinfection, bottles should be left to dry completely upside-down, in a bottle rack. Dry bottles should be closed with a lid and stored, protected from dust, until use.

# (v) Preparation of simple organophosphates, phosphonates and

#### thiophosphates:

#### **Preparation of organophosphates:**

In <u>organic chemistry</u>, organophosphates (also known as phosphate esters, or OPEs) are a class of <u>organophosphorus compounds</u> with the general structure  $O=P(OR)_3$ , a central phosphate molecule with <u>alkyl</u> or <u>aromatic substituents</u>.<sup>[11]</sup> They can be considered as <u>esters</u> of <u>phosphoric</u> <u>acid</u>. Organophosphates are best known for their use as <u>pesticides</u>.

Simple organophosphates can be synthesized via alcoholysis of phosphorus oxychloride (POCl<sub>3</sub>), esterification of phosphoric acid, or oxidation of phosphite esters. The most common industrial method is alcoholysis, where POCl3 reacts with alcohols to form organophosphates, often with the aid of catalysts like aluminum trichloride.

#### Alcoholysis of POCl<sub>3</sub>

Phosphorus oxychloride reacts readily with alcohols to give organophosphates. This is the dominant industrial route and is responsible for almost all organophosphate production. Aluminum trichloride or magnesium chloride are often employed as catalysts.

Reaction:  $POCl_3 + 3ROH \rightarrow O=P(OR)_3 + 3HCl$ 

Esterification of Phosphoric Acid:

**Reaction:**  $OP(OH)_3 + ROH \rightarrow OP(OH)_2(OR) + H_2O$ , followed by further esterification

Organophosphates are a class of compounds encompassing a number of distinct but closely related <u>function groups</u>. These are primarily the <u>esters</u> of <u>phosphoric acid</u> and can be mono-esters, di-esters or tri-esters depending on the number of attached <u>organic</u> groups (abbreviated as 'R' in the image below). In general man-made organophosphates are most often

triesters, while biological organophosphates are usually mono- or di-esters. The hydolysis of triesters can form diesters and monoesters.



In the context of pesticides, derivatives of organophosphates such as <u>organothiophosphates</u> (P=S) or <u>phosphorodiamidates</u> (P-N) are included as being organophosphates. The reason is that these compound are converted into organophosphates biologically.



In biology the esters of <u>diphosphoric acid</u> and <u>triphosphoric acid</u> are generally included as organophosphates. The reason is again a practical one, as many cellular processes involve the mono-, di- and tri- phosphates of the same compound. For instance, the phosphates of <u>adenosine (AMP, ADP, ATP)</u> play a key role in many metabolic processes.



#### Synthesis

#### **Alcoholysis of POCl3**

<u>Phosphorus oxychloride</u> reacts readily with <u>alcohols</u> to give organophosphates. This is the dominant industrial route and is responsible for almost all organophosphate production. <u>Aluminum trichloride</u> or <u>magnesium chloride</u> are often employed as <u>catalysts</u>.  $O=PCl_3 + 3 \text{ ROH} \rightarrow O=P(OR)_3 + 3 \text{ HCl}$ 

When aliphatic alcohols are used the HCl by-product can react with the phosphate esters to give <u>organochlorides</u> and a lower ester.

 $O=P(OR)_3 + HCl \rightarrow O=P(OR)_2OH + RCl$ 

This reaction is usually undesirable and is exacerbated by high reaction temperatures. It can be inhibited by the use of a base or the removal of HCl through <u>sparging</u>.

#### Esterification of phosphoric acid and P2O5

Esterifications of phosphoric acid with alcohols proceed less readily than the more common <u>carboxylic acid</u> esterifications, with the reactions rarely proceeding much further than the phosphate mono-ester. The reaction requires high temperatures, under which the phosphoric acid can dehydrate to form poly-phosphoric acids. These are exceedingly viscous and their linear polymeric structure renders them less reactive than phosphoric acid. Despite these limitations the reaction does see industrial use for the formation of monoalkyl phosphates, which are used as <u>surfactants</u>. A major appeal of this route is the low cost of phosphoric acid compared to phosphorus oxychloride.

 $OP(OH)_3 + ROH \rightarrow OP(OH)_2(OR) + H_2O$ 

 $\underline{P_2O_5}$  is the anhydride of phosphoric acid and acts similarly. The reaction yields equimolar amounts of di- and monoesters with no phosphoric acid. The process is mostly limited to primary alcohols, as secondary alcohols are prone to undesirable side reactions such as dehydration.

#### **Preparation of simple phosphonates:**

In <u>organic</u> chemistry, **phosphonates** or **phosphonic** acids are organophosphorus compounds containing C–PO(OR)<sub>2</sub> groups, where R is an organic group (<u>alkyl</u>, <u>aryl</u>). If R is hydrogen then the compound is a <u>dialkyl phosphite</u>, which is a different functional group. Phosphonic acids, typically handled as salts, are generally <u>nonvolatile</u> solids that are poorly <u>soluble</u> in <u>organic solvents</u>, but soluble in water and common <u>alcohols</u>.

Many commercially important compounds are phosphonates, including <u>glyphosate</u> (the active molecule of the herbicide <u>Roundup</u>), and <u>ethephon</u>, a widely used plant growth regulator. <u>Bisphosphonates</u> are popular drugs for treatment of <u>osteoporosis</u>.



Clodronic acid is a bisphosphonate used as drug treat osteoporosis. а to In biochemistry and medicinal chemistry, phosphonate groups are used as stable bioisosteres for phosphate, such as in the antiviral nucleotide analog, Tenofovir, one of the cornerstones of anti-HIV therapy. And there is an indication that phosphonate derivatives are "promising ligands for nuclear medicine.

#### Synthesis

Most processes begin with <u>phosphorous acid</u> (aka phosphonic acid, H<sub>3</sub>PO<sub>3</sub>), exploiting its reactive P–H bond.

#### From phosphonic acid

Phosphonic acid can be alkylated via the <u>Kabachnik–Fields reaction</u> or <u>Pudovik reaction</u> to give <u>aminophosphonate</u>, which are useful as <u>chelating</u> agents. One example is the industrial preparation of <u>nitrilotris(methylenephosphonic acid)</u>:

 $NH_3 + 3 H_3PO_3 + 3 CH_2O \rightarrow N(CH_2PO_3H_2)_3 + 3 H_2O$ 

Phosphonic acid also can be alkylated with <u>acrylic acid</u> derivatives to afford carboxyl functionalized phosphonic acids. This reaction is a variant of the <u>Michael addition</u>:

 $CH_2=CHCO_2R + 3 H_3PO_3 \rightarrow (HO)_2P(O)CH_2CH_2CO_2R$ 

In the <u>Hirao coupling</u> dialkyl <u>phosphites</u> (which can also be viewed as di-esters of phosphonic acid: (O=PH(OR)<sub>2</sub>) undergo a palladium-catalyzed coupling reaction with an <u>aryl halide</u> to form a phosphonate.

#### From phosphorus trichloride

<u>Vinylphosphonic acid</u> can be prepared by the reaction of PCl<sub>3</sub> and <u>acetaldehyde</u>:

 $PCl_3 + CH_3CHO \rightarrow CH_3CH(O^-)PCl^+_3$ 

This adduct reacts with acetic acid:

 $CH_3CH(O^-)PCl^+_3 + 2 CH_3CO_2H \rightarrow CH_3CH(Cl)PO(OH)_2 + 2 CH_3COCl$ 

This chloride undergoes dehydrochlorination to afford the target:

 $CH_3CH(Cl)PO(OH)_2 \rightarrow CH_2=CHPO(OH)_2 + HCl$ 

#### **Preparation of simple** Thiophosphates:

**Thiophosphates** (or phosphorothioates, PS) are <u>chemical compounds</u> and anions with the general <u>chemical formula</u>  $PS_{4-x}O^{3-}_{x}$  (x = 0, 1, 2, or 3) and related derivatives where organic groups are attached to one or more O or S. Thiophosphates feature tetrahedral phosphorus(V) centers.

Organothiophosphates are a subclass of <u>organophosphorus compounds</u> that are structurally related to the inorganic thiophosphates. Common members have formulas of the type  $(RO)_{3-x}(RS)_xPS$  and related compounds where RO is replaced by RS. Many of these compounds are used as <u>insecticides</u>, some have medical applications, and some have been used as <u>oil additives</u>.



Amifostine, which is used in cancer chemotherapy.



Chlorpyrifos, a popular insecticide.



Malathion, a popular insecticide.

#### **Synthesis**

Dithiophosphate has the formula  $[PO_2S_2]^{3-}$ , which has  $C_{2v}$  symmetry. Sodium dithiophosphate, which is colorless, is the major product from the reaction of <u>phosphorus</u> pentasulfide with <u>NaOH</u>:

 $P_2S_5 + 6 \text{ NaOH} \rightarrow 2 \text{ Na}_3PO_2S_2 + H_2S + 2 \text{ H}_2O$ 

**Dithiophosphoric acid** is obtained by treatment of barium dithiophosphate with <u>sulfuric acid</u>: Ba<sub>3</sub>(PO<sub>2</sub>S<sub>2</sub>)<sub>2</sub> + 3 H<sub>2</sub>SO<sub>4</sub>  $\rightarrow$  3 BaSO<sub>4</sub> + 2 H<sub>3</sub>PO<sub>2</sub>S<sub>2</sub>

Both  $Na_3PO_2S_2$  and especially  $H_3PO_2S_2$  are prone toward hydrolysis to their monothio derivatives.

# Part B: <u>Field Visit Report on the Preparation of Household Chemicals and its</u> <u>Applications</u>

#### **Introduction:**

Household chemicals are non-food chemicals that are commonly found and used in almost every place like homes, schools, offices and markets. They are designed particularly to assist cleaning house and back garden maintenance, pest control and for other hygiene purposes. Cleaning is the process of removing unwanted substances such as dirt, stains, infectious agents and other impurities from an object or environment So Cleaning is very essential to stay fit & healthy. Cleaning also protect the health of human beings by controlling the transmission of diseases.

For Cleaning we require different types of cleaning agents and disinfectants. A cleaning agent is not the instrument used for cleaning but rather the chemicals in liquid or solid forms actually does the cleaning. They are either natural or synthetic substances that are used to aid the cleaning process. The toxicity of a chemical substance depends on the strength and nature as well as how long the person is in contact with it. So it is important to know the various aspects of their uses. This course is designed to learn about Preparation of household cleaning agents and disinfectants.

#### **LEARNING OUTCOMES:**

- The students should learn fundamentals of household chemicals.
- ✤ The students should define various processes of house hold chemicals.
- The students should explain preparations, uses and reactions of household chemicals.
- The students should define preparation of different strength cleaning agents and hand Sanitizer.
- The students should prepare various composition of Liquid Detergent, liquid disinfectant and Dish washes Liquid Soap.
- The students should explain preparations and reactions of simple organophosphates, phosphonates and thiophosphates.
- ✤ The students should learn technology of detergent.
- Students should define herbal product, extraction techniques.
- Preparation of hand sanitizer of different flavor and fragrances.

MIDNAPORE CITY COLLEGE

#### **Objectives:**

- > To study the Preparation of Household Chemicals and its Applications.
- > To collect information regarding various use and its chemical hazards.
- > To identify traditional use in public health and safety.

#### **Observation:**

Name of place: Insert image

#### **Description:**

**Chemical constituents:** 

**Traditional use:** 

# Paper Title: Inorganic Chemistry –II & Organic Chemistry-II (Practical) Paper Code: CEMMIN03P

# **Course Content:**

#### **Inorganic Chemistry Laboratory**

1. Estimation of Fe(II) and Fe(III) in a given mixture.

2. Estimation of Fe(III) and Cu(II) in a mixture.

#### **Organic Chemistry Lab**

A. Separation, based upon solubility, by using common laboratory reagents like water (cold, hot), dil. HCl, dil. NaOH, dil. NaHCO3, etc., of components of a binary solid mixture; purification of any one of the separated components by crystallization.

The composition of the mixture may be of the following types: Benzoic acid/p-Toluidine; p-Nitrobenzoic acid/p-Aminobenzoic acid; p-Nitrotolune/p-Anisidine; etc.

B. Identification of a pure organic compound

Solid compounds: oxalic acid, tartaric acid, succinic acid, resorcinol, urea, glucose, benzoic acid and salicylic acid.

Liquid Compounds: methyl alcohol, ethyl alcohol, acetone, aniline, dimethylaniline, benzaldehyde, chloroform and nitrobenzene.

### **Organic Chemistry Lab**

#### **1. Separation**

**Basic Concept:** Synthesis of organic compounds often involves reactions that result in the formation of a mixture of products. When some parts of a tree are extracted in an attempt to isolate a targeted compound, invariably a mixture of compounds appear. Isolation of the desired compound, therefore, needs exclusion of the unwanted companions. This constitutes the principle of separation of a mixture of compounds into its components. Wide variety of techniques, including solvent extraction, fractional crystallisation and steam distillation, chromatography etc are employed for this purpose. An increase in the number of components in the mixture requires a technique with enhanced sophistication and clinical accuracy for separation. In the following section separation of a mixture of only two components utilising the difference in solubility of components in a suitable solvent is discussed.

If one component dissolves in a particular solvent and the other does not, the mixture may be shaken with the appropriate solvent to dissolve the soluble component. Filtration separates the insoluble component from the soluble one. Then the soluble component has to be regenerated from its solution. Acidification, basification, evaporation etc may be needed for this regeneration purpose.

#### **General Procedure of Separation:**

 Test the solubility of the supplied Binary Mixture by adding a pinch of the sample in (-2 mL) each of the following two solvents and shake well to find the solvent in which a part of the mixture is soluble.

Solvents: (i) 10% NaOH solution, (ii) 6 (N) HCI.

- 2. Add (-25 mL) of the appropriate solvent (identified in step 1) to the rest of the Mixture taken in a 100 mL beaker and stir well with a glass rod for 15 20 minutes.
- 3. Filter under suction and collect the filtrate in a 250 mL beaker.
- 4. Wash the residue obtained in step 3 with the same solvent used in step 2. Then wash the residue with water, dry and collect as the 'Insoluble Part'

5. Precipitate the soluble component from the filtrate obtained in step3 either by acidifying with conc. HCl (if 10% NaOH solution is used for dissolution) or basifying with 20% NaOH solution (if 6(N) HCl solution is used for dissolution) as the case may be. Filter off the precipitate, wash with water, dry and collect the solid as the soluble parts.

#### **Reporting of Separation of Components of Unknown Binary Solid Mixture:**

Chemical required	Apparatus required
6(N) HCl	Test tube
Conc. HCl	100 ml beaker
10% NaOH	Measuring
	cylinder

Chemical	Apparatus
required	required
20% NaOH	Glass rod
Distilled water	Suction filtering
	apparatus with
	small Buchner
	funnel

**TYPE I [Mixture of Acidic (PART A) and Neutral (PART B) Components] :** 0.2 g of the given binary mixture (preferably in powdered form, if crystalline sample is given, it should be ground thoroughly with the pestle to make it powder) is vigorously shaken separately with 2-3 mL of each of 6(N) HCl and 10% NaOH solution at room temperature for 5 minutes.

6(N) HCl	10% NaOH
Insoluble	Partly Soluble



Melting point of part A= Melting point of part B=

**TYPE II [Mixture of Basic (PART A) and Neutral (PART B) Components] :** 0.2 g of the given binary mixture (preferably in powdered form, if crystalline sample is given, it should be ground thoroughly with the pestle to make it powder) is vigorously shaken separately with 2-3 mL of each of 6(N) HCl and 10% NaOH solution at room temperature for 5 minutes.

10% NaOH	6(N) HCl
Insoluble	Partly Soluble



Melting point of part A=

Melting point of part B=

#### **Identification of a Pure Organic Compound**

The sample supplied may be any one of the following types:

- 1. Pure solids
- 2. Pure liquid

The following preliminary tests will be performed with the given sample;

# Solid compounds

#### **Oxalic acid**

#### Physical characteristics and preliminary test :

1. State	: Solid
2. Colour	: Colorless
3. Texture	: Crystaline

MIDNAPORE CITY COLLEGE

3. Odour	: Odourless
4. Solubility	: Soluble in cold water
5. Litmus	: Blue litmus turns red
6. Action of heat	: Volatilises completely without
	charring when heated slowly and
	a sublime is formed but
	decomposed with evolution of
	gases when heated strongly
7. Ignition test	: Blue non-sooty flame

# <u>Succinic acia</u>

- 1. State
- 2. Colour
- 3. Texture
- 3. Odour
- 4. Solubility

- : Solid
- : Colorless
- : Crystaline
- : Odourless
- : Soluble in cold water

5. Litmus	: Blue litmus turns red
6. Action of heat	: Melts and then boils giving off
	extremely irritating vapors. A
	sublimate is forme d at the
	cooler part of the test tube
7. Ignition test	: Blue non-sooty flame

Experiment	Observation
<b>1.</b> H <sub>2</sub> SO <sub>4</sub> test: Few crystals of O.S. is warmed conc. $H_2SO_4$ and then heated saongly.	1. O.S dissolve without charring. Slight chairring occurs on strong heating and the solution turn brown with evolution of SO <sub>2</sub> .
<ol> <li>Calcium chloride test: An equal volume of Calcium chloride solution is added to neutral solution of O.S. The mixture is shaken vigorously inner side of the test tube is scratched with glass rod and finally boiled.</li> <li>Silver nitrate test : AgNO<sub>3</sub> solution is added to neutral solution of O.S.</li> <li>Denige's test (C.T.): Denige's reagent is added to neutral solution or aqueous solution of O.S. and then a drop of 2% KMnO4 solution is added.</li> </ol>	<ol> <li>A white ppt. soluble in acetic acld.</li> <li>Buff-colored ppt soluble in dil. HCl.</li> <li>No ppt. in the cold but the color of permanganate presists.</li> </ol>

# **Resorcinol**

1. State	: Solid
2. Colour	: Colorless
3. Texture	: Crystaline
4. Odour	: Odourless
5. Solubility	: Soluble in cold water
6. Litmus	: feebly acidic
7. Action of heat	: Melts and volatiles
8. Ignition test	: Yellow sooty flame

Experiment	Observation
<b>1.</b> $H_2SO_4$ test: Few crystals of O.S. is warmed with conc. $H_2SO_4$ .	1. O.S dissolve without charring.
<b>2. Ferric chloride test</b> : A few drops of FeCl <sub>3</sub> solution is added to aqueous solution of O.S.	2. Bluish-violet color.

<b>3. Bromine writer test:</b> Bromine writer is	3. White ppt.
added to aq. solution of O.S. and shaken.	4. Black ppt of silver mirror.
<b>4. Tollen's test:</b> Aqueous solution of O.S is added	
to Tollen's reagent and warmed.	5 An intense green fluorescein
5. Fluorescein test: A small quantity of O.S. is	5. 7 m mense green nuoreseem.
mixed with succinic acid and the test performed.	

# <u>Urea</u>

I hysical characteristics and premining y test	Physical	characteristics	and	preliminary	test	:
--	----------	-----------------	-----	-------------	------	---

1. State	: Solid
2. Colour	: Colorless
3. Texture	: Crystaline
4. Odour	: Odourless
5. Solubility	: Soluble in cold water
6. Litmus	: Neutral
7. Action of heat	: Melts and gives off characteristic
	smell of NH <sub>3</sub>
8. Ignition test	: Blue non-sooty flame

Experiment	Observation
1. Nitrous acid test: Dil. HCl and dil. NaNO2	1. Effervescence with evolution of
solution is added to aq. Solution of O.s in a test	gas.
tube and shaken.	
	2. A white crystalline ppt.
<b>2.</b> Nitric acid test: 2 mL of conc. HNO3 is added	
to aqueous concentrate solution of O.S.	3. A white crystalline ppt.
<b>3. Oxalic acid test:</b> A concentrate solution of	4. A pink or violet colour develops.
oxalic acid is added to conc. solution of O.S.	r r
<b>4. Biuret test:</b> About 0.5 g of O.S. is heated gently	
to melt. The heating is continued till the molten	
mass solidifies at once when test tube is taken out	
of flame. If there is no solidification the process	
is repeated. The residue is cooled and dissolved in	
I ml of NaOH solution by warming and shaking.	
The solution is cooled again and a drop or two of	
a very dilute copper sulphate solution is added to	
it.	

## **Benzoic acid**

- 1. State
- 2. Colour
- 3. Texture

- : Solid
- : White
- : Plate-shaped or needle shaped

	Crystal
4. Odour	: Odourless
5. Solubility	: Insoluble in cold water but
	Soluble in hot water
6. Litmus	: Blue litmus turns red
7. Action of heat	: Melts and sublimates
8. Ignition test	: Yellow sooty flame

Experiment	Observation
1. H <sub>2</sub> SO <sub>4</sub> test: A small amount of O.S. is	1. Dissolves without charring.
warmed with conc. $H_2SO_4$ .	
2. Soda-lime test: A small amount of	2. Characteristics smell of benzene.
O.S. is heated in a hard glass test tube	
with sodalime.	3. White ppt. dissolves on boiling but
3. Denige's test: Denige's reagent is	appears on cooling.
added to neutral solution of O.S.	4. (a) a buff-colored ppt.
<b>4. Ferric chloride test:</b> (a) A drop or two	(b) buff-colored ppt dissolves with
of freshly prepared solution of Fecl <sub>3</sub> is	appearance of white ppt.
added to neutral solution of O.S. (b)	·FF ··································
Dilute HCl is added to it.	5 characteristics fruity small of ester
<b>5. Esterification test:</b> About 0.1 g of O.S	5. characteristics fruity shield of ester.
and 2 mL of dehydrated alcohol are taken	
in a clean test tube and shaken well. Then	
a few drops of conc. H <sub>2</sub> SO <sub>4</sub> is added and	
the test tube is waremd in a hot water	
bath. The test tube is cooled and content	
of the test tube is poured into 50 mL of	
dilute sodium cxarbonate solution taken	
in a beaker.	

# Salicylic acid

Physical characteristics and preliminary test :

- 1. State
- 2. Colour
- 3. Texture
- 4. Odour
- 5. Solubility
- 6. Litmus
- 7. Action of heat
- 8. Ignition test

- : Solid
- : White
- : needle-shaped Crystal
- : Odourless
- : Sparingly soluble in cold water but readily Soluble in hot water
- : Blue litmus turns red
- : Melts and sublimes
- : Yellow sooty flame

MIDNAPORE CITY COLLEGE

Experiment	Observation
<b>1. H<sub>2</sub>SO<sub>4</sub> test:</b> A small amount of O.S. is warmed with conc. H <sub>2</sub> SO <sub>4</sub> .	1. The solid dissolved, charring occurs after some time, solution darkens and gas evolve.
<ol> <li>Soda-lime test: A small amount of O.S. is heated in a hard glass test tube with sodalime.</li> <li>Denige's test: Denige's reagent is added to neutral solution of O.S.</li> <li>Ferric chloride test: A drop or two of freshly prepared solution of Fecl<sub>3</sub> is added to neutral solution of O.S. (a) Dilute HCl is added to it. (b) acetic acid is added to another part.</li> <li>Oil of winter test: About 0.1 g of O.S and 1 mL of methanol and few drops of conc. H<sub>2</sub>SO<sub>4</sub> are warmed in aclean dry test tube.</li> </ol>	<ol> <li>Characteristics smell of phenol.</li> <li>White ppt. dissolves on boiling but appears on cooling.</li> <li>An intense violet color         <ul> <li>(a) The color discharges</li> <li>(b) The color persists</li> </ul> </li> <li>A characteristics pungent fragnant odor which intensifies on pouring the mixture to dil. Sodium carbonate solution.</li> </ol>

# Tartaric acid

- 1. State
- 2. Colour
- 3. Texture
- 4. Odour
- 5. Solubility
- 6. Litmus
- 7. Action of heat
- 8. Ignition test

- : Solid
- : Colorless
- : Crystalline
- : Odourless
- : Soluble in cold water
- : Blue litmus turns red
- : Charrs and gives off smell of burnt sugar
- : Blue non-sooty flame

Experiment	Observation
1. H <sub>2</sub> SO <sub>4</sub> test: A small amount of O.S. is	1. charrs immediately.
warmed with conc. $H_2SO_4$ .	
2. Calcium chloride test: Excess of	2. A white crystalline ppt. soluble in hot
calcium chloride solution is added to	dilute acetic acid.

neutral solution of O.S. shaken and	
inside of the test tube is scratched with a	
glass rod.	$A$ (a) $\Delta$ white crystalline ppt (b)
4. Silver nitrate test : (a) AgNO <sub>3</sub>	shining silver mirror forms at the
solution is added to neutral solution of	inner of the test tube
O.S. (b) Addition of AgNO3 solution is	
continued till the precipitation is	
complete. Then dil. NH4OH is added	
dropwise with constant shaking till the	
ppt. almost and not completely dissolved	5. Permanganate color discharges
as revealed by slight turbidity.	immediately without turbidity.
5. Denige's test (C.T.): Denige's reagent	
is added to neutral solution or aqueous	
solution of O.S.	

# **Glucose**

- 1. State
- 2. Colour
- 3. Texture
- 4. Odour
- 5. Solubility
- 6. Litmus
- 7. Action of heat
- 8. Ignition test

- : Solid
- : Colorless
- : Crystalline
- : Odourless
- : Soluble in cold water
- : Neutral
- : Charrs and gives off smell of burnt sugar
- : Blue non-sooty flame

Experiment	Observation
1. H <sub>2</sub> SO <sub>4</sub> test: A small amount of O.S. is	1. No charring in cold but darkens on
warmed with conc. $H_2SO_4$ .	heating immediately.
2. Sodium hydroxide test: 2 ml of 3%	
solution of NaOH is added to conc.	2. The solution turns first yellow and
Solution of O.S. and the mixture is	then raddish brown.
heated. The mixture is then acidified with	
dil. HNO <sub>3</sub> .	
3. Tollen's test : 5 ml of aqueous	3. Black or grev precipitate or shining
solution of o.s is added to aeual volume	silver mirror.
of Tollen's reagent taken in a clean test	
tube and then the test tube is placed in a	1 Brick red precipitate
beaker of boiling water for few minutes.	4. Drick red precipitate.

4. Fehling's test: The sample is added to	
a mixture of Fehling A and Felling B and	
heated for 1 minute.	5. The white ppt. turns salmon-pink in
5. Lead acetate test: Lead acetate	color.
solution is added to solution of o.s and	
the mixture is boiled for few seconds.	
Then dil. NH <sub>4</sub> OH solution is added	
dropwise till just sufficient to produce a	
permanat white ppt. then the mixture is	
boiled again.	

# Liquid compounds Methyl Alcohol

CH<sub>3</sub>-OH

# Physical characteristics and preliminary test :

1. State	: Liquid
2. Colour	: Colorless
3. Odour	: Pungent but faint alcoholic
4. Miscibility	: Miscible With water
5. Litmus	: Neutral
6. Action of heat	: Volatilises
7. Ignition test	: Blue non-sooty flame

Experiment	Observation
<b>1. Oxidation test :</b> A Cu-spiral is made	1.
repeatedly red hot and introduced into 2	
ml of O.S. kept in a test tube dipped in a	
beaker containing cold Water. The	
solution is then divided into four pan ts.	
(a) Schiffis test : One part is added to	(a) Pink color which deepens slowly.
Schiff s.	
(b) Tollen's test: Another part is	(b) Bright mirror of silver on the side of
added to Tollen's reagent and heated	the test tube.
in a water-bath.	
(d) <b>Resorcinol test:</b> A mixture of 0.5 ml	(c) Intense indigo blue color
of oxidised liquid and a drop of 0.5% aq.	(c) mense mengo blue color.
Solution of resorcinol is added carefully	
down the side of the test tube held in an	
inclined position containing 2ml of conc.	
$H_2SO_4.$	
<b>2. Denige's reagent:</b> To about 5ml of aq.	
Solution of o.s. in a porcelain basin	

placed on ice-water, 2-3 ml of 2.5% soln. of KMnO <sub>4</sub> soln. is added followed by few drops of conc. $H_2SO_4$ when a brown color developes. The brown color is destroyed by addition of saturated soln. of oxalic acid. A freshly prepared soln.of schiffs reagent is added dropwise to it with stirring.	<ul><li>(d) A reddish -violet ring appears at the junction of two liquid layers. A white ppt. changing to reddish-violet appear in the aq. Layer above the ring after a little while.</li><li>2. A violet color develops deepening on standing.</li></ul>
<b>3.</b> Oil of wintergreen test: 0.5 ml of o.s. and 3 drops of conc. H2SO4 are added to 0.5 g of salicylic acid and heated for one minute and finally the reaction mixture is poured into 50 ml of water taken in a	
beaker.	3. Characteristics smell of methyl salicylate.

# Ethanol CH<sub>3</sub>CH<sub>2</sub>OH

Physical characteristics and preliminary test :

1. State	: Liquid
2. Colour	: Colorless
3. Odour	: Pungent but faint alcoholic
4. Miscibility	: Miscible With water
5. Litmus	: Neutral
6. Action of heat	: Volatilises
7. Ignition test	: Blue non-sooty flame

Experiment	Observation
<b>1. Oxidation test :</b> A Cu-spiral is made	1.
repeatedly red hot and introduced into 2	
ml of O.S. kept in a test tube dipped in a	
beaker containing cold Water. The	
solution is then divided into four pan ts.	
(a) Schiffis test : One part is added to	(a) Pink color which deepens slowly.
Schiff s.	
(b) Tollen's test: Another part is	(b) A silver mirror is disposited on the
added to Tollen's reagent and heated	walls of the test tube.
in a water-bath.	
(c) To another part equal volume of very	

dilute solution of sodium nitroprusside solution is added followed by a few drops of NaOH solution.	(c) Wine-red color develops.
(d) To another part a very dilute solution of sodium nitroprusside solution is added followed by a few drops of piperidine.	(d) A deep red color develops.
<b>2. Ethyl ester formate:</b> About 0.2 ml of o.s an equal volume of glacial acetic acid or 0.2g of fused sodium acetate and a little conc. H2SO4 are taken in a test tube and warmed in a water-bath for about 5min. the mixture is then poured into water taken in a backer to which a little	2. Pleasant fruity smell.
<ul> <li>3. Iodoform test: (a) To 2 ml of aq. Soln. of o.s. an equal volume of a conc. Solution of iodine in potassium iodide is added and then NaOH solution is added dropwise with stirring until violet color</li> </ul>	3. (a) Characteristics smell of methyl salicylate.
of iodine dis appears. Then the mixture is warmed and cooled under tap with shaking. (b) The iodoform test is repeated by adding NH4OH instead of NaOH.	(b) No ppt. of iodoform forms.

# Chloroform CHCl<sub>3</sub>

### Physical characteristics and preliminary test :

- 1. State
- 2. Colour
- 3. Odour
- 4. Miscibility
- 5. Litmus
- 6. Action of heat
- 7. Ignition test

- : Liquid
- : Colorless
- : Sweet
- : Immiscible With water
- : Neutral
- : Volatilises
- : Yellow sooty flame

Experiment	Observation
------------	-------------

<b>1. Hydrolysis test:</b> o.s. is treated with aq.	1.Curdy white ppt. soluble in NH4OH
KOH or NaOH solution and boiled. The	but reappears on addition of conc.
mixture is acidified with conc. HNO3	HNO3.
and then AgNO3 is added.	
2. Fehlings test: o.s and Fehlings	2. A yellowish red ppt.
solution is warmed with constant	
shaking.	3. A brilliant reddish colouration is
3. Resorcinol test: A little powdered	developed in aq. Laver.
resorcinol and few drops of O.S. is taken	
in a basin, about 1 ml of conc. solution of	
NaOH is added to it, then the mixture is	
warmed gently.	
	4. Intolerable obnoxious smell of
4. Carbylamine test: Few drops of O.S.	carbylamine.
aniline and alc. KOH are warmed in a diy	
test tube and the ensuing gas is smelt by	
placing the thumb at the mouth of the test	
tube and then holding tfie thumb near the	
nostrils.	

**Nitrobenzene** 



Physical characteristics and preliminary test :

1. 9	State
------	-------

- 2. Colour
- 3. Odour
- 4. Miscibility
- 5. Litmus
- 6. Action of heat
- 7. Ignition test

- : Liquid
- : Pale yellow
- : Characteristic smell of bitter almonds
- : Immiscible With water
- : Neutral
- : Volatilises without leaving any residue
- : Yellow sooty flame

Experiment	Observation
------------	-------------

1 Caustic soda test · About 0.2 ml of	1 The color darkens
1. Caustic soua test . About 0.2 III of	1.The color darkens.
O.S. is heated with conc. solution of	
caustic soda.	
2. Reduction and diazocoupling test:	
0.5 ml is reduced by tin or zinc and dil.	
HCl for 5 minutes. The decanted	2. Brilliant scarlet-red dye.
solution is cooled and very dilute cold	
solution of NaNO2 is added to it. The	
diazotised soln. is added to cold alkaline	
soln. of. β-naphtol.	3. Grev or black ppt
3. Muliken-Barker test: An aq.	5. Grey of black ppt
Ethanolic solution of few drops of given	
sample is boiled with a pinch of zinc dust	
and little solid NH4Cl. The mixture is	
cooled and filtered into Tollen's reagent.	

# Acetone CH<sub>3</sub>COCH<sub>3</sub>

Physical characteristics and preliminary test :

- 1. State
- 2. Colour
- 3. Odour
- 4. Miscibility
- 5. Litmus
- 6. Action of heat
- 7. Ignition test

- : Liquid
- : Colorless
- : Pleasant ethereal odour
- : Miscible With water
- : Neutral
- : Completely Volatilises without leaving any residue
- : Blue non-sooty flame

Experiment	Observation
<b>1. Oxidation test :</b> A small amount of O.S. is added to 2-4-dinitrophenylhydrazine solution.	1. Yellow ppt. on mixing.
<b>2. Denige's test:</b> An equal volume of Denige's reagent is added to an aqueous solution of O.S and the test tube is then kept in a boiling water bath for few minutes.	2. Heavy white ppt. of double compound of acetone and basic mercuric sulphate.
<b>3. Legal's test:</b> A few drops of a very dil. solution of sodium nitroprusside is added to 2-3 ml of aqueous solution of O.S. followed	3. A ruby red (reddish-purple) color develops which disappears on warming but reappears on cooling.
<ul><li>by a few drops of dil. NaOH solution.</li><li>4. Iodoform test: (a) iodoform test is performed with I2 in KI solution and NaOH</li></ul>	<ul> <li>4. (a) Yellow ppt. of CHI<sub>3</sub> having characteristics sweeet smell.</li> <li>(b) Yellow ppt. of CHI<sub>3</sub> having</li> </ul>

solution. (b) the iodoform test is repeated by	characteristics sweeet smell.
adding NH <sub>4</sub> OH solution instead of NaOH	
solution.	



# Physical characteristics and preliminary test :

1. State	: Liquid
2. Colour	: Colorless when freshly distilled
	but turns brown on exposure to
	light and air
3. Odour	: Characteristic aromatic smell
4. Miscibility	: Immiscible with water but soluble
	dil.HCI
5. Litmus	: Neutral
6. Action of heat	: Completely Volatilises without
	leaving any residue
7. Ignition test	: Yellow sooty flame

Experiment	Observation
1. Bleaching powder test: A few drops of	1. Purple-violet colouration.
bleaching powder solution is added to a dilute	
solution of O S. in HCl.	
2. Potassium dichromate test: A drop of	2. Intense blue colour develops.
O.S. is added to 5 to 6 drops of conc. $H_2SO_4$	1
taken in a spot plate and the mixture is stirred	
well with glass rod. Then a drop of K <sub>2</sub> Cr <sub>2</sub> O <sub>7</sub>	2 Prilliant scorlat rad dya
solution is added to it.	5. Diffiant scallet-fed uye.
<b>3. Diazocoupling test:</b> 5-6 drops of O.S are	
dissolved in dil. HCl in a test tube and cooled	
in ice-water. Then 3-4 drops of very dilute ice	
cold solution of sodium nitrite are added to it.	
Then the solution is added to 1mL of ice-cold	
alkaline solution of $\beta$ -naphtol.	

# Dimethylanilinebenzaldehyde



Physical characteristics and preliminary test :

1. State	: Liquid
2. Colour	: Colorless when freshly distilled
	but turns brown on exposure to
	light and air
3. Odour	: Characteristic bad smell
4. Miscibility	: Immiscible with water but soluble
	dil.HCI
5. Litmus	: Neutral
6. Action of heat	: Completely Volatilises without
	leaving any residue
7. Ignition test	: Yellow sooty flame

# **Special Tests**

Experiment	Observation
1. Bleaching powder test: A few drops of	1. Purple-violet colouration.
bleaching powder solution is added to a dilute	
solution of O S. in HCl.	
2. Potassium dichromate test: A drop of	2. Intense blue colour develops.
O.S. is added to 5 to 6 drops of conc. $H_2SO_4$	1
taken in a spot plate and the mixture is stirred	
well with glass rod. Then a drop of K <sub>2</sub> Cr <sub>2</sub> O <sub>7</sub>	2 Prilliont georlat rad due
solution is added to it.	5. Brinant scarlet-red uye.
3. Diazocoupling test: 5-6 drops of O.S are	
dissolved in dil. HCl in a test tube and cooled	
in ice-water. Then 3-4 drops of very dilute ice	
cold solution of sodium nitrite are added to it.	
Then the solution is added to 1mL of ice-cold	
alkaline solution of $\beta$ -naphtol.	

# Benzaldehyde CH<sub>3</sub>CHO

1. State	: Liquid
2. Colour	: Colorless
3. Odour	: Characteristic smell of bitter almonds
4. Miscibility	: Immiscible with water
5. Litmus	: Neutral
6. Action of heat	: Volatilises
7. Ignition test	: Yellow sooty flame

# **Special Tests**

Experiment	Observation
<b>1. 2,4-D.N.P test:</b> A few drops of 2,4-D.N.P	1. reddish yellow ppt. forms immediately
solution is added to 1 mL alcoholic solution	simply on mixing.
given sample.	
bleaching powder solution is added to a dilute	
solution of O S. in HCl.	2. Black ppt.
2. Tollen's test: An alcoholic solution of the	11
given sample is added to 1 mL of Tollen's	
reagent warmed gently.	2 A intense green soleration
3. Malachite green test: 0.5 ml of O.S is	5. A intense green coloration.
heated with 1ml of dimethyl aniline and a	
small bit of anhydrous ZnCl2 in a dry test	
tube for one minute. The leucobase produced	
is oxidised with lead dioxide in a solution of	
acetic acid and excess conc. HCl is added.	

#### **Inorganic Chemistry Laboratory**

# **Experiment 1: Estimation of Fe(II) and Fe(III) in a given mixture.**

AIM:

To estimate the amount of ferrous and ferric ions in a given mixture by making use potassium dichromate

Chemicals: potassium dichromate, Stannous chloride, Mercurous chloride, phenyl amine

**Preparation of solutions**:

Standard potassium dichromate:

Weigh out accurately weighted 0.245 gm potassium dichromate into a 100 ml standard flask, dissolve in little amount of water, and make up the solution with distilled water. Shake the flask well for uniform concentration.

#### Stannous chloride solution:

Dissolve 150 g of SnCl2.2H2O in 100 mL concentrated HCl and dilute with water to 1 L

#### Mercuric chloride solution:

Dissolve 60 g of HgCl2 in hot water and dilute to 1 litre.

#### Acid mixture:

Dissolve 150 ml of concentrated H2SO4 and 50 ml of H3PO4 in a reagent bottle and stopper it.

#### **Phenylamine:**

Dissolve 1 gm of di Phenyl amine in 100ml of concentrated H2SO4.

#### **Procedure:**

#### **Standardisation of total iron:**

Make up the given mixture of  $Fe^{+2}$  and  $Fe^{+3}$  given in 100 ml standard flask up to the mark withdistilled water and shake the flask well for uniform concentration. Pippette out 20 ml of the iron solution in to a conical flask, Add 100 ml concentrated HCl and heat the solution to bolingtill the colour of the solution changes yellow. Add stannous chloride drop by drop to the hot solution till the yellow colour disappears. Cool the solution under tap and add 100 ml mercuricchloride in one portion silky white ppt of mercurous is obtained. Add 5 ml of acid mixture and 3- 4 drops of diphenylamine indicator .Titrate the solution against standard potassium dichromate taken in a burette till the green colour changes to blue violet. Repeat the titration toget concurrent values. Let the titre value be x ml.

#### Standardization of ferrous iron:

Pipette out 20 ml of the total iron solution into a 250 ml conical flask, add 5 ml of acid mixture and 2 drops of di-phenylamine indicator. Titrate the resulting solution with potassium dichromate taken in a burette till a blue colour is obtained. Note the burette readings and repeat the titration for concurrent values. Let the titre value be y ml

#### **CALCULATIONS:**

Normality of potassium dichromate N1= Weight of potassium dichromate  $\times 10/49$ Normality of total iron = V1N1=V2N2 V<sub>1</sub>= Volume of dichromate  $V_2 = Volume of dichromate consumed for$ total ironN2 = Normality of total iron N1 = Normality of dichromate N2 = V1N1/V2= titre value be x ml  $\times$ N1/20 Amount of total iron = N2  $\times$ 55.85 = p gms /100 ml 10 Normality of ferrous iron  $Fe^{+2}=$  $V_1N_1 = V_3N_3V_1 = V_0$  of dichromate N1 = Normality ofdichromate V3 = Volumeof  $Fe^{+2}$  solution N3 = Normality of Fe<sup>+2</sup> solution N3 = V1N1 = titre value be y ml×N1/20V 3 Amount of Ferrous iron = N3  $\times$ 55.85= q gms /100ml 10

#### **RESULT:**

Amount of Ferrous iron =Amount of total iron-amount of  $Fe^{+2} = p - q/100$  gms

THE C

#### Experiment 2: Estimation of Fe(III) and Cu(II) in a mixture

Principle : For the estimation of individual components in Fe(III) + Cu(II) ions mixture, both dichromatometry as well as iodometry are employed. Since both Fe(III) and Cu(II) ions can oxidize iodide to iodine, for preferential oxidation of iodide by Cu(II) ions only, Fe(III) ions are at first to be masked using ammonium bifluoride. This reduces the formal potential of Fe(III)/Fe(II) couple lower than the  $E^0$  value = 0.54V (of  $I_2/I^-$  couple) and thus Fe(III) ions can no longer oxidize I- to I2. As Cu(II) ions don't form stable complex with fluoride ions but with iodide ions (following HSAB principle), Cu(II) ions now oxidize iodide to iodine and itself get reduced to Cu(I) to form white precipitate of CuI. The precipitation effect increases the formal potential of Cu(II)/Cu(I) couple higher than  $E^0$  value = 0.54V of  $I_2/I^-$  couple and thus increases the spontaneity of the reaction. The stoichiometric quantity of iodine thus liberated is then titrated with standard thiosulphate using starch indicator. Near the end point the deep blue colour of starch-iodine complex is faded, 1% NH4SCN is added, shaken well as the blue colour of the solution deepens due to libration of adsorbed iodine on CuI which is desorbed on transformation of CuI to CuSCN. Thus, a creamy white suspension results at the end point of titration when blue colour of starch-iodine complex is completely discharged.

$$2Cu^{2+} + 4I^{-} \rightleftharpoons 2CuI \text{ (white ppt)} + I_2$$

$$I_2 + 2S_2O_3^{2-} \rightleftharpoons 2I^{-} + S_4O_6^{2-}$$

$$2CuI + 2SCN^{-} \rightarrow 2CuSCN + I_2$$

$$2 \text{ mole of } Cu^{2+} \equiv 1 \text{ mole of } I_2 \equiv 2 \text{ mole of } Na_2S_2O_3$$

$$1 \text{ mole of } Na_2S_2O_3 \equiv 1/2 \text{ mole of } I_2 \equiv 1 \text{ gm-atom of } Cu^{2+}$$

$$1000 \text{ mL 1(N) } Na_2S_2O_3 \equiv 63.54 \text{ g } Cu^{2+}$$

Fe (III) ions are estimated by separating them from Cu(II) ions through precipitation of Fe(OH)3 using ammonium hydroxide in presence of ammonium chloride and then re-dissolving it in 6N HCl followed by Al foil or SnCl2 reduction method to form Fe(II) ions which is then titrated with standard K2Cr2O7, in presence of phosphoric acid or ammonium bifluoride masking agent using BaDS indicator. At the end point, the green solution changes to red-violet. Al-foil reduction :  $Fe^{3+} + Al \rightarrow 3Fe^{2+} + Al^{3+}$ 

SnCl<sub>2</sub>-HCl reduction :

$$2Fe^{3+} + Sn^{2+} \rightarrow 2Fe^{2+} + Sn^{4+}$$

 $SnCl_2 + HgCl_2 \rightarrow SnCl_4 + 2Hg_2Cl_2$  (silky white precipitate)  $SnCl_2 + Hg_2Cl_2 \rightarrow SnCl_4 + 2Hg$  (colloidal mercury)

**Overall reaction :** 

$$Cr_2O_7^{2-}+14H^++6Fe^{2+} \rightleftharpoons 2Cr^{3+}+7H_2O+6Fe^{3+}$$

$$Cr_2O_7^{2-} \equiv 6Fe^{2+}, 1/6 Cr_2O_7^{2-} \equiv 1 Fe^{2+} \equiv 1 eqv \equiv 1 Fe^{3+}$$

1000 mL 1(N) Cr<sub>2</sub>O<sub>7</sub><sup>2-</sup> = 55.85g Fe<sup>3+</sup>

#### **Chemicals required:**

- () Standard (N/20) K2Cr2O7 solution and all the second states and the state
- (1) 4N H2SO4 solution
- (3) Syrupy phosphoric acid or ammonium bifluoride solid
- ④ 0.1% BaDS indicator
- (5) Concentrated HCl
- (6) Al foil or SnCl<sub>2</sub> solution [150g of SnCl<sub>2</sub>·2l·l<sub>2</sub>O in 500mL conc. HCL and dilute to 1 L with distilled water].
- (7) 10% HgCl<sub>2</sub> solution
- (8) 20 mL 10% KI solution
- (9) 1% starch solution
- (10) NHACI
- (ff) 1:1 (v/v) Ammonia
- (1) Fe(III) + Cu(II) stock solution: (i) Dissolve 24g ferric alum in 500 mL 2(N) sulphuric acid (ii) Dissolve 12g of CuSO4 · 5H2O salt solution in in 500mL 2(N) sulphuric acid. Add (i) and (ii) mixtures to prepare stock solution
- (B) 1% NH<sub>4</sub>SCN (w/v) in distilled water.
- (N/20) Na<sub>2</sub>S<sub>2</sub>O<sub>3</sub> solution (See section 2.5 for details)
- (15) Concentrated Nitric acid

Procedure :

Step-1 : Preparation of 250 mL standard N/20 K2Cr2O7 solution.

(Follow the procedure as discussed in section- 2.4.1)

Step-2 : Standardization of Na2S2O3 solution by standard (N/20) K2Cr2O7 solution.

(Follow the procedure as discussed in section- 2.5.1, step-2)

Step-3 : Estimation of Cu(II) ions : 25 mL of supplied stock solution is pipetted out in a 250 mL conical flask and neutralized by 1:1 NH4OH till a faint permanent bluish white turbidity appears. Now add 2 g of NH4HF2 to dissolve the precipitate, followed by 2 g of solid KI or 20 mL 10% KI solution and dilute to 100 mL. The solution turns brown immediately due to iodine liberation which is then titrated with standard Na2S2O3 solution to straw yellow colour. Now add 2 mL of 1% starch solution, as the solution turns deep blue. Continue titration with standard Na2S2O3 solution till the blue colour of starch is faded away. Further add 10 mL 1% NH4SCN, shake well as the blue colour of the solution deepens and continue the titration till blue colour is completely discharged and creamy white solution results. Record the titre value and repeat the experiment.

Step-4 : Estimation of Fe(III) ions : 25 mL of Fe(III) ions stock solution is pipetted out in a 500 mL beaker and diluted to 100 mL by water. Now 2 g of NH4Cl is added, heated to boiling and 1:1 (v/v) ammonia is added dropwise till smell of ammonia persists. Reddish brown precipitate to ferric hydroxide appears and is allowed to settle down on a hot asbestos board, till the supernatant liquid shows no further precipitation of ferric hydroxide. Filter hot the resultant solution through Whatman-41 filter paper, wash with 1% NH4Cl aqueous solution containing few drops of ammonia. Now the precipitate is dissolved in minimum volume of hot, 1:1 (v/v) HCl, followed by hot water and is collected in the same beaker where initial

precipitation was carried out. (Re-precipitation process is carried out to remove unwanted, adsorbed Cu(II) ions).

Reduce the solution by Al foil and heat the solution carefully (avoid boiling), shake the solution continuously till the yellow colour is discharged.

[The reduction can also be carried out by adding the SnCl<sub>2</sub> reagent dropwise with constant stirring. When the yellow colour of the solution is fully discharged, add 1 drop of extra diluted SnCl<sub>2</sub> (1 : 1 HCl solution) and cool the solution under tap to room temperature. Then add 10 mL 10% HgCl<sub>2</sub> to the solution all at a time and swirl the solution and wait for 10 minutes. A silky white precipitate appears due to formation of Hg<sub>2</sub>Cl<sub>2</sub>. Excess reagent may otherwise transform Hg<sub>2</sub>Cl<sub>2</sub> to colloidal mercury that may interfere in the titration and impart error to analysis.]

After reduction to Fe(II) ions, cool the solution and add 2g  $NH_4HF_2$  or 5 mL phosphoric acid followed by 4-5 drops BaDS indicator. Titrate with standard (N/20)  $K_2Cr_2O_7$  solution till the colour of the solution just changes from green to violet. Record the titre value.

#### **Results** :

Estimation of Cu(II) ions :

1. Preparation of 250 mL standard (N/20) K2Cr2O7 solution :

Weight of K <sub>2</sub> Cr <sub>2</sub> O <sub>7</sub>	Weight of K <sub>2</sub> Cr <sub>2</sub> O <sub>7</sub> to be	Strength of K <sub>2</sub> Cr <sub>2</sub> O <sub>7</sub>
taken (g)	taken (g)	solution
w <sub>k</sub>	0.6129	$\mathbf{F} = \left(\frac{w_k}{0.6129 \times 20}\right) (\mathbf{N})$

2. Standardization of Na<sub>2</sub>S<sub>2</sub>O<sub>3</sub> solution by standard (N/20) K<sub>2</sub>Cr<sub>2</sub>O<sub>7</sub> :

No. of titrations	Volume of F(N) standard K <sub>2</sub> Cr <sub>2</sub> O <sub>7</sub> solution	Volume of thiosulphate S <sub>2</sub> O <sub>3</sub> <sup>2-</sup> consumed	Mean volume of consumed thiosulphate solution (mL)
1	25.0 mL	<u>X</u> mL ·	$-V_2 = (X + Y)/2$
.2	25.0 mL	YmL	

N.B.: If consecutive readings do not differ by ± 0.1 mL, then a third set of titration to be performed.

So, strength of  $Na_2S_2O_3$  solution (S<sub>t</sub>) =  $[25 \times F / V_2]$  (N)

3. Estimation of Cu2+ ions :

No. of titrations	Volume of stock solution	Volume of thiosulphate consumed	Mean volume of consumer thiosulphate solution (mL)
1	25.0 mL	A mL	$V_t = (A + B)/2$
2 , 1 .	25.0 mL	BmL	

N.B.: If consecutive readings do not differ by ± 0.1 mL, then a third set of titration to be performed.

■ Calculation : 1000 mL 1 (N) Na2S2O3 = 63.54 g Cu

 $V_t \text{ mL } S_t(N) \text{ Na}_2 S_2 O_3 \text{ solution} == [(63.54/1000) \times V_t \times S_t] \text{ g of } Cu$ Now, 25 mL stock solution contains ==  $[(63.54/1000) \times V_t \times S_t]$  g of Cu Therefore, 1000 mL - - - =  $[(40 \times 63.54)/1000 \times V_t \times S_t] = h \text{ g of } Cu$  Conclusion : The given/supplied solution contains  $\dots h \dots g/L$  of  $Cu^{2+}$  ion.

No. of	Volume of stock	Volume of standard [F(N)]	Mean volume of standard
titrations	solution	K <sub>2</sub> Cr <sub>2</sub> O <sub>7</sub> solution.	[F(N)] K <sub>2</sub> Cr <sub>2</sub> O <sub>7</sub> solution.
1	25 mL	V <sub>Fe</sub> mL	V <sub>Fe</sub> mL

■ Calculation: 1000 mL 1 (N) K<sub>2</sub>Cr<sub>2</sub>O<sub>7</sub> = 55.85 g Fe

$$V_{\rm Fe} \, \text{mL F(N)} \, \text{K}_2 \text{Cr}_2 \text{O}_7 \equiv [55.85 \times V_{\rm Fe} \times \text{F}/1000] \, \text{g Fe} = z \, \text{g of Fe}$$

Now , 25 mL stock solution contains z g of Fe

So 1000 mL stock solution contains (40 x z) g of Fe

■ Conclusion : The supplied solution contains ......g/lit of Fe<sup>3+</sup> ion

\*\*\*\*\*\*